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82319 Starnberg • GERMANY

info@pari.de • www.pari.com

# Reprocessing instructions: All valid procedures

## Including optional procedures

Model: PARI PEP® S System (Type 018)

Model: PARI LC PLUS (Type 022)

Model: PARI LC SPRINT® (Type 023)

Model: PARI Mask soft (Type 041)

Model: PARI BABY mask (Type 041)

Model: PARI SMARTMASK® (Type 041/078)

Model: PARI BOY® (Type 130)

Model: PARI COMPACT2 (Type 152)

## **Read the document and instructions for use**

Always observe the instructions for use supplied with the PARI product. All necessary information for pretreatment and preparation are described in the instructions for use. Read the instructions carefully before using the product. Follow all instructions and safety directions. Keep the instructions in a safe place. The document contains additional information with which pretreatment and preconditioning can also be achieved.

## **Document validity**

PARI PEP® S System (Type 018)

PARI LC PLUS® (Type 022)

PARI LC SPRINT® (Type 023)

PARI Maske soft (Type 041)

PARI BABY Maske (Type 041)

PARI SMARTMASK® (Type 041/078)

PARI BOY® (Type 130)

PARI COMPACT2 (Type 152)

## **Request the latest version**

You can order this document from PARI GmbH.

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Moosstrasse 3, 82319 Starnberg, Germany. Technical and design changes and printing errors reserved. The same applies to illustrations.

## **Disclaimer**

This document describes the reprocessing of PARI products and optional accessories. Therefore, it is possible that not all components and sections for reprocessing described will be relevant to your PARI product. This document also describes and illustrates features not present in your PARI product because they are, for instance, country-specific and/or optional. When using the systems, products and functions, the applicable country-specific regulations must be observed.

## **Trade marks**

Registered trade marks of PARI GmbH Spezialisten für effektive Inhalation in Germany and/or other countries:

BOY®, LC SPRINT®, PARI®, PEP®

## **Trade marks**

Registered trade marks of PARI Pharma GmbH in Germany and/or other countries:

PARI SMARTMASK®

# TABLE OF CONTENTS

<b>IMPORTANT INFORMATION</b> .....	4
<b>REPROCESSING</b> .....	7
Reprocessing cycles.....	8
Limits of reprocessing.....	10
Preparation for reprocessing .....	11
Reprocessing of parts.....	11
Reprocessing the compressor.....	18
Reprocessing the connection tubing .....	20

# 1 IMPORTANT INFORMATION

## 1.1 Safety and warning instructions

This document contains important information, safety instructions and precautionary measures. The user must follow these in order to guarantee safe operation of this PARI product.

Only use this PARI product as described in this document and the instructions for use of the corresponding product.

The instructions for use of the inhalation solution used must also be followed.

### Labelling and classification of warning instructions

Safety-critical warnings are categorised in this document according to the following hazard levels:

#### **DANGER**

DANGER indicates a hazardous situation which will lead to very severe injuries or death if it is not avoided.

#### **WARNING**

WARNING indicates a hazardous situation which can lead to very severe injuries or death if it is not avoided.

#### **CAUTION**

CAUTION indicates a hazardous situation which can lead to mild or moderate injuries if it is not avoided.

#### **NOTE**

NOTE indicates a hazardous situation which can lead to material damage if it is not avoided.

## Life-threatening situation from electrocution

**DANGER**

### Life-threatening situation from electrocution

The compressor is an electrical device powered by mains voltage. It has been designed so that no live parts are accessible. However, in unfavourable ambient conditions, or if the compressor or power cord is damaged, this protection may no longer be provided. There may then be a risk of contact with live parts. This in turn may lead to an electric shock.

• To avoid this danger, comply with the following instructions:

- Every time before use, make sure that the compressor housing, the power cord and the power plug are undamaged. The compressor must not be operated
  - if the housing, the power cord or the power plug is damaged,
  - if a fault is suspected following a fall or similar.
- Never leave the compressor unattended while it is in use.
- Plug the compressor into an easily accessible socket. It should be possible to quickly unplug the power plug at all times.
- Switch the compressor off and disconnect the power plug from the mains socket immediately:
  - if it is suspected that the compressor or the power cord might have been damaged (e.g. after the compressor has fallen, or if there is a smell of burning plastic)
  - if a malfunction occurs during operation
  - before cleaning and maintaining the device
  - immediately after use
- Keep the power cord away from domestic animals (e.g. rodents). They may damage the insulation on the power cord.

## Danger of a device defect

To avoid a device defect, comply with the following instructions:

- Make sure that the local supply voltage matches the voltage data marked on the compressor identification label.
- To disconnect the compressor from the socket, always pull the power cord out by the power plug, not the cable.
- Make sure that the power cord is never kinked, pinched or jammed. Do not pull the power cord over sharp edges.
- Keep the compressor and the power cord away from hot surfaces (e.g. stove top, electric fire, open fire). Direct heat may damage the compressor housing or the insulation on the power cord.

## Hazard due to small parts which can be swallowed

The product contains small parts. Small parts can block the airways and lead to a choking hazard. Keep all components of the product out of the reach of babies and infants at all times.

## Impairment of treatment due to electromagnetic interference

Use only original spare parts and original accessories from PARI. The use of third-party products can result in increased emissions of electromagnetic interference or reduced interference resistance of the PARI compressor.

## Impairment of the therapy by disregarding the reprocessing instructions

### CAUTION

#### Impairment of the therapy

Exceeding the reprocessing instructions can result in damage to the product. This can impair the therapy.

Falling short of the instructions can result in insufficient reprocessing. This can increase the risk of infection.

- Comply with the instructions about the reprocessing limits, temperature, holding time, and concentration of the chemical used.

## Reporting serious incidents

Report serious adverse incidents to the manufacturer and to the competent authority.

## 2 REPROCESSING

### CAUTION

#### **Risk of infection due to cross-contamination in the case of a change in patients**

If a product is used for more than one patient, there is a risk that germs may be transmitted from one patient to the next.

- Clean, disinfect and sterilise all separated parts before every patient change.
- Replace the connection tubing or carry out mechanical cleaning and disinfection of the connection tubing [see: Reprocessing the connection tubing, page 20].

### CAUTION

#### **Risk of infection due to growth of bacteria**

Inadequate disinfection encourages the growth of bacteria and thus increases the risk of infection.

- During disinfection, comply with the specified holding time and the concentration of the chemical used.
- When using a disinfectant, ensure that the device is clean and operating properly. Allow the disinfection to continue until the disinfectant switches off automatically or the minimum disinfection time stated in the instructions for use of the disinfectant has elapsed. Do not switch the device off prematurely.

### CAUTION

#### **Risk of infection by residual germs**

If there is dirt on the parts, germs capable of reproduction may remain despite the sterilisation process. As a result, there is a danger of infection.

- Clean, disinfect, and dry all parts thoroughly before sterilising.
- Use only validated procedures for cleaning and disinfection.

### CAUTION

#### **Risk of infection due to moisture**

Moisture encourages the growth of bacteria.

- Ensure that the components are allowed to dry sufficiently after each processing step.

## NOTE

### Risk of damage to plastic parts!

Plastic melts if it comes into contact with hot surfaces.

- During thermal disinfection, you must ensure that the water level in the container is sufficient, in order that the individual parts do not come into contact with hot surfaces.

*Always use drinking water for reprocessing in a home environment.*

 *In professional health care institutions, drinking water is sufficient for precleaning. For all other reprocessing steps and processes, use deionised water with a low microbiological load (at least drinking water quality).*

 *If the recommended application period is exceeded significantly, the plastic parts may take on the smell of the medium used.*

Please observe the instructions for use for the chemicals used.

Wash your hands thoroughly before every reprocessing.

## 2.1 Reprocessing cycles

Perform the reprocessing steps as per the specified reprocessing cycles.

### Reprocessing cycles in a home environment

Nebuliser components	– Clean immediately after every use – Disinfect once per week <sup>1</sup>
Mask components	– Clean immediately after every use – Disinfect once per week <sup>1</sup>
Elastic band	Cleaning visible dirt
PARI PEP S parts	– Clean immediately after every use – Disinfect once a week <sup>1</sup>
Compressor housing	– Cleaning in the event of visible soiling and before each patient change – Disinfect by wiping before each patient change
Air filter	Replace after 200 operating hours (approx. 1 year) <sup>2</sup>

1) For the therapy of patients at risk, disinfect the individual parts once per day. Take note of the additional information for patients at risk [see: Reprocessing for patients at risk, page 10].

2) Air filters are included in every PARI Year Pack.

# Reprocessing cycles in professional healthcare institution

## **SINGLE PATIENT USE**

Nebuliser components	– Clean immediately after every use – Disinfect once per week <sup>1</sup>
Mask components	– Clean immediately after every use – Disinfect once per week <sup>1</sup>
Elastic band	Cleaning visible dirt
PARI PEP S parts	– Clean immediately after every use – Disinfect once a week <sup>1</sup>
Compressor housing	Cleaning visible dirt
Air filter	Replace after 200 operating hours (approx. 1 year) <sup>3</sup>

## **BEFORE A CHANGE OF PATIENTS**

Nebuliser components	– Cleaning – Disinfection – Sterilisation
Mask components	– Cleaning – Disinfection – Sterilisation
Connection tubing	Mechanical cleaning with disinfection
Elastic band	Replacing the elastic band
PARI PEP S parts	– Cleaning – Disinfection – Sterilisation <sup>4</sup>
Compressor housing	– Cleaning – Wipe disinfection
Air filter	Replace after 200 operating hours (approx. 1 year) <sup>3</sup>

3) Air filters are included in every PARI Year Pack.

4) The nose clip cannot be sterilised.

## Reprocessing for patients at risk

### CAUTION

#### Danger of infection for patients at risk

For patients at risk, airway infections represent an increased risk of deterioration of general health, because these patients are particularly vulnerable to residual germs. Patients at risk include cystic fibrosis patients, patients with immunosuppression or immunodeficiency, and vulnerable patient groups.

- If you are a patient at risk, disinfect the individual parts once per day during the therapy.
- If you are unsure whether you are a patient at risk, consult with specialist medical personnel before use.

## 2.2 Limits of reprocessing


### Limits of reprocessing in a home environment

Nebuliser and accessories, disinfection	300 processing cycles, max. 1 year
Mask, disinfection	300 processing cycles, max. 1 year
PARI PEP S System, disinfection	300 processing cycles, max. 1 year

### Limits of reprocessing in a professional healthcare institution

Nebuliser and accessories, disinfection	300 reprocessing cycles, max. 1 year
Nebuliser and accessories, sterilisation	100 processing cycles, max. 1 year
Connection tubing	50 reprocessing cycles, max. 1 year
Mask, disinfection	300 processing cycles, max. 1 year
Mask, sterilisation	100 processing cycles, max. 1 year
PARI PEP S System, disinfection	300 processing cycles, max. 1 year
PARI PEP S System, sterilisation <sup>4</sup>	100 processing cycles, max. 1 year

## 2.3 Preparation for reprocessing


Procedure	
Compressor	<ul style="list-style-type: none"><li>• Detach the connection tubing from the compressor.</li></ul>
Nebuliser	<ul style="list-style-type: none"><li>• Detach the connection tubing from the nebuliser.</li><li>• Detach the mouthpiece from the nebuliser.</li><li>• Make sure that all residual volume is removed from the nebuliser.</li><li>• Dismantle the nebuliser into its individual parts.</li><li>• Carefully pull the blue exhalation valve out of the slot in the mouthpiece. The exhalation valve must still be attached to the mouthpiece.</li></ul> 
Mask	<ul style="list-style-type: none"><li>• Disconnect all mask components from the nebuliser.</li><li>• If necessary, detach the bend piece (45° angle) from the mask.</li><li>• Detach the PARI BABY bend from the mask if applicable.</li><li>• If necessary, detach the elastic band from the mask.</li></ul>
PEP S	<ul style="list-style-type: none"><li>• If necessary, disconnect the PARI PEP S system from the nebuliser.</li><li>• Dismantle the PARI PEP S system into its individual parts.</li><li>• Open the cap on the PEP S connecting element.</li></ul>


## 2.4 Reprocessing of parts

The following products can be cleaned, disinfected and sterilised according to the procedure described below:

- PARI nebuliser and PARI accessories
- PARI mask soft
- PARI BABY mask
- PARI SMARTMASK
- PARI BABY bend
- PARI bend piece (45° angle)
- PARI PEP S System

The connection tubing and elastic band of the mask must be treated separately.

## Procedure

### Step 1: Precleaning

All individual parts must be precleaned immediately after use.

EQUIPMENT:

– Drinking water temperature of about 15 °C

PROCEDURE:

- Rinse all parts used for 2 minutes in running drinking water.

### Option A - Manual cleaning

#### EQUIPMENT:

- Drinking water temperature of about 40 °C
- Standard commercial washing-up liquid<sup>5</sup>
- Container having sufficient capacity

#### PROCEDURE:

- Unless otherwise specified by the manufacturer of the washing-up liquid, add about 1 teaspoon of washing-up liquid to 3 l warm drinking water.
- Place all the parts in the washing-up water.  
Application time: 5 minutes
- Occasionally move the parts back and forth.
- In case of visible soiling, use a medium-soft brush (e.g. a toothbrush) which is used exclusively for this purpose.

#### RINSING:

- Rinse off all parts thoroughly in running drinking water at about 15 °C for 3 minutes.

### Option B - Machine cleaning

The individual parts can be cleaned in a standard household dishwasher provided it is connected to a mains water supply of drinking water quality.

#### PROCEDURE:

- Place all components in the crockery basket so that no water can collect in them.
- Select a program with at least 50 °C.

 *Do not clean the individual components together with very dirty dishes.*

### Cleaning the elastic band

- Clean the elastic band as necessary with warm drinking water and a little dishwashing liquid.

 *The elastic band cannot be disinfected or sterilised.*

## Procedure

Home environment	Step 3: Disinfection	<p><b>Option A - Thermal disinfection using boiling water</b></p> <p>EQUIPMENT:</p> <ul style="list-style-type: none"> <li>– Clean cooking pot</li> <li>– Drinking water</li> </ul> <p>PROCEDURE:</p> <ul style="list-style-type: none"> <li>• Place all the individual parts in water at a rolling boil for at least 5 minutes.</li> </ul> <hr/> <p><b>Option B - Using a standard thermal disinfectant for baby bottles</b></p> <p>EQUIPMENT:</p> <ul style="list-style-type: none"> <li>– Thermal disinfectant with a runtime of at least 6 minutes.</li> </ul> <p>PROCEDURE:</p> <p>Regarding the performance of the disinfection, the duration of the disinfection procedure and the quantity of water required, follow the instructions for use of the disinfectant you are using.</p>
In professional health care institutions	Step 2: Cleaning and disinfection	<p><b>Mechanical cleaning and disinfection</b></p> <p>EQUIPMENT:</p> <ul style="list-style-type: none"> <li>– Neodisher® MediClean forte (concentration: 0.5%)</li> <li>– Deionised water</li> <li>– Cleaning and disinfection device as per ISO 15883.</li> </ul> <p>PROCEDURE:</p> <ul style="list-style-type: none"> <li>• Programme for cleaning and thermal disinfection (at least A0 = 3000) as per manufacturer's instructions.</li> </ul>
Drying		<p>Ensure that there is no residual moisture remaining in the components. If necessary:</p> <ul style="list-style-type: none"> <li>• Shake the water out of all of the parts.</li> <li>• Place all parts on a dry, clean and absorbent surface and allow them to dry completely.</li> <li>• Using a standard thermal disinfectant for baby bottles: Leave all individual parts in the closed thermal disinfectant for max. 24 hours until the next use.</li> </ul>

## Procedure

<p style="writing-mode: vertical-rl; transform: rotate(180deg);">In professional health care institutions</p>	<p style="writing-mode: vertical-rl; transform: rotate(180deg);">Step 3: Sterilisation</p>	<p>For sterilisation, use a process compliant with the ISO 17665 series of standards.</p> <p>PROCEDURE:</p> <p>Sterilisation temperature and holding time:</p> <ul style="list-style-type: none"> <li>– Temperature: 134 °C</li> <li>– Holding time: minimum 3 minutes to maximum 5 minutes</li> </ul>	
		<p><b>Information on sterilising a PARI mask soft:</b></p> <p>Always use the corresponding mask stabiliser when sterilising this mask type, because otherwise the mask may lose its shape under the effects of high temperatures.</p>	
		<ul style="list-style-type: none"> <li>• Insert the mask stabiliser in the mask as shown in the illustration.</li> <li>• Pack all the disassembled parts in a sterile barrier system according to ISO 11607-1 (e.g. foil-paper packaging).</li> <li>• Carry out the sterilisation in a steam steriliser in accordance with the manufacturer's instructions.</li> </ul>	
<p style="writing-mode: vertical-rl; transform: rotate(180deg);">Visual inspection</p>	<p>Inspect all product components after each cleaning, disinfection or, where applicable, sterilisation. Replace defective, deformed or seriously discoloured parts.</p>		
<p style="writing-mode: vertical-rl; transform: rotate(180deg);">Storage</p>	<p>Store this product as described below:</p> <ul style="list-style-type: none"> <li>– Dry and dust-free, e.g. in a clean, lint-free cloth (e.g. a tea towel)</li> <li>– Protected against contamination, where necessary (e.g. using optional sterile packaging)</li> </ul>		

## 2.5 Reprocessing the compressor

### DANGER

#### Life-threatening situation from electrocution

Liquids can conduct electricity, thereby posing a risk of electric shock.

- Before starting to clean the compressor, always switch it off, and disconnect the power plug from the mains socket.

#### NOTE

#### Danger of device fault due to liquid penetration


If liquids get into the interior of the compressor, this may cause a fault in the device.

- Never immerse the compressor in water.
- Never clean the compressor in running water.
- Never spray any liquids onto the compressor or the power cord.
- If liquid gets into the compressor, it must not be used under any circumstances. Before starting the compressor again, contact the manufacturer or distributor.

### Procedure

#### Step 1: Cleaning

- Wipe the outer surface of the housing with a clean, damp cloth.

#### Step 2: Disinfection

For disinfection, use a standard, alcohol-based disinfectant (e.g. isopropanol). When applying and dosing the disinfectant, it is essential to follow the instructions for use for the product.

- If the compressor is visibly dirty, clean the compressor before disinfecting it.
- Moisten a cloth with the disinfectant.
- Wipe off the outer surface of the housing thoroughly with the cloth.
- Let the disinfectant dry completely on the surfaces.

The procedure was validated in Europe using a disinfectant suitable for use on plastics as per the DGHM and/or the VAH list: Propanol/Isopropanol  
Validated using Incidin® liquid.

## Procedure

### Visual inspection


- Inspect all product components after each cleaning and disinfection.
- Replace defective, deformed or seriously discoloured parts.

The air filter must be checked at regular intervals (after every 10th use). If it is discoloured brown or grey, or if it is damp or clogged, it must be replaced.

**The air filter cannot be cleaned and then reused!**

#### Removing the air filter

- Pull the filter holder out of the compressor. Use a small screwdriver, for example, to prise the filter holder carefully out of the compressor.


#### Replacing the air filter

- **NOTE!** For operating your compressor, use only air filters provided by the manufacturer or distributor. If air filters not designed for the compressor are used, it may be damaged.  
Pull the old air filter out of the filter holder, and fit the new filter in its place.
- Insert the filter holder back in the compressor.

### Storage

- **⚠ CAUTION! Electrical devices connected to the power supply represent a potential source of danger.**  
When storing the compressor, always pull the power plug out of the socket.

- **NOTE!** Do not wind the power cord around the compressor. If the power cord is wound or bent very tightly, the wires inside the cord may break. The power cord will then be unusable.  
Wind the power cord up loosely.
- Secure the cable end in the cable holder.
- Wrap all individual components in a clean, lint-free cloth (e.g. a tea towel).
- Store the product in a dry, dust-free place.

## 2.6 Reprocessing the connection tubing

Dry the connection tubing after each use.

The maximum operating life of the connection tubing is 1 year.

Procedure	
In professional health care institutions (when changing between patients)	<b>Mechanical cleaning and disinfection</b> <b>EQUIPMENT:</b> The method has been validated in Europe using: <ul style="list-style-type: none"><li>– Neodisher® MediClean forte (concentration: 0.5%)</li><li>– Deionised water</li><li>– Cleaning and disinfection device as per ISO 15883.</li><li>– Special baskets for Miele instrument dishwasher</li><li>– Compressed air source for drying</li></ul> <b>PROCEDURE:</b> <ul style="list-style-type: none"><li>• Programme for cleaning and thermal disinfection (at least A0 = 3000) as per manufacturer's instructions.</li></ul>
<b>Drying</b>	<ul style="list-style-type: none"><li>• Connect the connection tubing to a compressed air source (compressor or central medical gas supply system).</li><li>• Switch the compressed air source on.</li><li>• Leave the compressed air source running until all the moisture in the tube has been removed.</li></ul>
<b>Visual inspection</b>	<ul style="list-style-type: none"><li>• Inspect all product components after each cleaning and disinfection.</li><li>• Replace defective, deformed or seriously discoloured parts.</li></ul>
<b>Storage</b>	<ul style="list-style-type: none"><li>• Store the product as described below:<ul style="list-style-type: none"><li>– Dry and dust-free, e.g. in a clean, lint-free cloth (e.g. a tea towel)</li><li>– Protected against contamination, where necessary (e.g. using optional sterile packaging)</li></ul></li></ul>

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## **Copyright**