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Instruction for Use

INQUAir Compressor

Registered and Production address: 2nd floor, Block A, No.58, Fumin South Road, Dalang Town, Dongguan City, Guangdong Province, China
Date of production: see label
After-sale service: INQUA GmbH Germany
Email: inquair@inqua.de

Registered and manufactured by: Dongguan City Tianchang Electro-Mechanical Manufacturing Co., Ltd

1 of 14

Index

1. Name.....	3
2. Scope.....	3
3. Basic Structure.....	3
4. Technical Specifications.....	3
5. Contraindications.....	4
6. Packaging.....	4
7. Storage and Shipping Conditions.....	4
8. Accessories.....	4
9. How to Use.....	5
10. Cautions.....	5
11. Marking on Outer Package.....	6
12. Electromagnetic Compatibility.....	6
13. Maintenance.....	11
14. Troubleshooting.....	12
15. Warranty Card.....	13
16. Packing List (in each box).....	14

We reserve the right to change the appearance and/or the technology of this product.

Please keep this instruction for future reference.

- See page 7 for the electromagnetic compatibility of this product
- Do not drop, detach or modify this product. Avoid direct sunlight, humid and dusty environment.

11. Marking on Outer Package

Sign				
Notes	This Side Up	Handle with care	Max Stacks	Avoid Rain

Marking on Device

Sign			
Notes	Electric Shock Proof - Type 2	Electric Shock Proof Intensity - BF Application	Water Proof Level - IPXO

12. Electromagnetic Compatibility

- INQUAir Compressor model ANBB26024 is consistent with the requirements of the YY0505-2012 standard electromagnetic compatibility.
- Users should install and use this product according to the information on the electromagnetic compatibility provided in the box.
- Portable and mobile radio frequency communication device may affect the performance of the INQUAir Compressor model ANBB26024. Please avoid strong electromagnetic interference such as cellphone and microwave oven during use.
- See appendix 1 for Guidelines and manufacturer's Declaration

6 of 14

The recommended safe distance between portable and mobile radio frequency communication device and the compressor model ANBB26024

Compressor model ANBB26024 is expected to be used in an electromagnetic environment where the radio frequency transmission interference is controlled. According to the max. rated power output of the communication device, buyer or user can maintain a min. distance between portable and mobile radio frequency transmitter and compressor model ANBB26024.

Max. Rated Power Output in watt	Safe distance in different frequency in meter			
	150 kHz ~ 80 MHz	80 MHz ~ 800 MHz	800 MHz ~ 2.5 GHz	
0.01	0.12	0.12	0.23	
0.1	0.38	0.38	0.73	
1	1.2	1.2	2.3	
10	3.8	3.8	7.3	
100	12	12	23	

For the max. rated power output not listed above, it is recommended the distance to be d in meters determined by the formula in corresponding frequency column. P is the max. rated power output by transmitter manufacture. The unit is watt.

Note 1: at 80MHz and 800MHz frequency, higher frequency range formula is adopted.

Note 2: the guidelines may not be applicable for certain situations. Electromagnetic transmission is affected by the absorption and reflection of buildings, objects and human body.

13. Maintenance

a. daily maintenance

Wipe clean the plastic housing with clean wet cloth. Place the INQUAir Compressor in its box when not in use.

b. replace the filter cotton

Filter cotton is to ensure the air flowing in is clean and uncontaminated. It is recommended to replace the filter cotton with a new piece when its color starts to change. There is a piece of spare filter cotton in the box. More spare filter cotton can be purchased from your supplier.

11 of 14

1. Name

Product Name: INQUAir Compressor
Model Number: ANBB26024
Dimensions: 168x120x211mm
Weight: 1400g

2. Scope

This product works with single patient use nebulizer kits including air tube, nebulizer cup, mask or mouthpiece and provides compressed air to aerosolize liquid medication for the use of patient.

3. Basic Structure

INQUAir Compressor mainly consists of a motor compressor and nebulizer kit PARI LC Sprint SPA or PARI LC Plus.

4. Technical Specifications

Electric Shock Proof Type	Type II	Input	220V/50Hz
Electric Shock Proof Intensity	BF Application	Power	180VA
Water Proof Level	IPX0	Fuse	250V/2A
Thermal Cut-Out Temperature	135°C		

Working flow rate	>4L/min	Working pressure	>82.7kpa
Maximum pressure	207kpa~380kpa	Atomization rate	≥0.20g/min
Noise level	≤65dB(A Weighted)	Residue	≤1.0g
Particle Concentration	Diameter < 5µm Particle Concentration ≥70	Continuous working hours	Under normal temperature, powered by AC power source, work continuously for 4 hours

Voltage Fluctuation/Flicker Emission GB 17625.2	Consistent	
---	------------	--

Guidelines and Manufacturer Declaration – Electromagnetic Interference Resistance


Compressor model ANBB26024 is expected to be used under the following electromagnetic environment.

Interference Resistance Test	IEC 60601 Test Electrical Level	Consistent Electrical Level	Electromagnetic Environment - Guideline
Static Current Discharge GB/T 17626.2	±6 kV Contact Discharge ±8 kV Air Discharge	±6 kV Contact Discharge ±8 kV Air Discharge	The surface should be wooden, concrete or tile. If the surface is covered by synthetic material, the relative humidity should be at least 30%. The power source should be of the quality that is acceptable by typical commercial or hospital environment
Electric Fast Transient GB/T 17626.4	±2kV to power cable ±1kV to input/output cable	±2kV to power cable	The power source should be of the quality that is acceptable by typical commercial or hospital environment
Surge GB/T 17626.5	±1 kV Cable to cable ±2 kV Cable to ground	±1 kV Cable to cable	The power source should be of the quality that is acceptable by typical commercial or hospital environment

7 of 14

c. how to replace the filter

- Open the cap
- Take out the filter cotton piece
- Put on a new piece and close the cap

14. Troubleshooting

No.	Issue	Possible Cause	Troubleshooting
1	Abnormal noise	Check whether the filter cotton is installed	Install the filter cotton by instructions
2	Low pressure and flow rate	Blocked filter cotton Life span exceeded	Reinstall filter cotton by instructions Replace the product
3	Machine won't start	Fuse blown Switch connector loose	Please contact us Please contact us
4	Machine makes noise but there is no atomization	Clam nylon screw loose	Please contact us
5	Water drop inside air tube	Air is too humid	Connect the air tube to INQUAir Compressor. Press the air outlet of the tube with your finger then release. Repeat the process.

Reminder:

If the INQUAir Compressor needs to be detached or its parts need replacing due to issues listed above, please contact your reseller for repair, electric schematics, part list or technical drawings.

Contact information:

INQUA GmbH
GERMANY
inquair@inqua.de

15. Warranty Card

After Sales Service

Product Name : INQUAir Compressor
Model : _____
Factory Serial Number : _____
Date of Purchase : _____

Notes:

Under normal circumstance and when in proper use, this product is covered by one-year limited warranty starting from the date of purchase, during which for any quality issue that is caused by non-artificial factor,

Transportation and storage	Atmospheric pressure: 500hPa~1060hPa	Working environment	Atmospheric pressure: 700hPa~1060hPa
	Atmospheric temperature: -20°~55° Atmospheric humidity: 10%-93% RH		Atmospheric temperature: 5°~40° Atmospheric humidity: 30%-80%RH

5. Contraindications

Patients who are not suitable for aerosol therapy are not allowed to use this product.

6. Packaging

Color box with one INQUAir Compressor unit, one nebulizer kit for single patient use, one IFU, one nebulizer kit IFU and one bag of filter cotton inside.

7. Storage and Shipping Conditions

Do not stack heavy object on it. Avoid direct sunlight and rain water. Do not put it in cold storage. Keep it indoor, at room temperature and in dry environment.

8. Accessories

Nebulizer kit (see the IFU for nebulizer).

Interference Resistance Test	IEC 60601 Test Electrical Level	Consistent Electrical Level	Electromagnetic Environment - Guideline
The temporary decrease, cutoff and change of the voltage on the power input cable GB/T 17626.11	<5 % U _i , lasting for 0.5 period (at U _i , >95% temporary decrease) 40 % U _i , lasting for 5 periods (at U _i , 60% temporary decrease) 70 % U _i , lasting for 25 periods (at U _i , 30% temporary decrease) <5 % U _i , lasting for 5s (at U _i , >95% temporary decrease)	<5 % U _i , lasting for 0.5 period (at U _i , >95% temporary decrease) 40 % U _i , lasting for 5 periods (at U _i , 60% temporary decrease) 70 % U _i , lasting for 25 periods (at U _i , 30% temporary decrease) <5 % U _i , lasting for 5s (at U _i , >95% temporary decrease)	Grid power source should be of the quality that is acceptable by typical commercial or hospital environment. If compressor model ANBB26024 is still to function during a power shortage, it is recommended an uninterruptible power supply or battery is used to drive the compressor. The power frequency magnetic field should be of the quality that is acceptable by typical commercial or hospital environment
Power Frequency Magnetic Field (50Hz) GB/T 17626.8	3A/m	3A/m, 50Hz	The power frequency magnetic field should be of the quality that is acceptable by typical commercial or hospital environment

8 of 14

Guidelines and Manufacturer Declaration – Electromagnetic Interference Resistance

user is entitled to return, repair or replace their products for free. User can contact their reseller with the invoice or proof of purchase for the service.

The following situations are not covered by warranty:

1. Damaged or deformed by force
2. Faults caused by water leakage or exposure to rain water
3. Water or medicine is sucked into the pump by user so that the INQUAir Compressor doesn't work properly
4. Wearables such as filter cotton and fuse are not covered


16. Packing List (in each box)

Name	Quantity
INQUAir Compressor	1
Mask	1
Mouthpiece	1
Nebulizer Cup	1
Air Tube	1
Product Qualification Certificate	1
Instruction for Use	2
Filter Cotton	1

Manufacturer: Dongguan City Tianchang Electro-Mechanical Manufacturing Co., Ltd
Address: 2nd floor, Block A, No.58, Fumin South Road, Dalang Town, Dongguan City, Guangdong Province, China

9. How to Use

- Place the INQUAir Compressor on a firm and flat surface.
- Connect the INQUAir Compressor to power then switch it on. Let it run unloaded for one minute to warm up, then switch it off (ignore this step if the environment temperature exceeds 5 degrees Celsius).
- Connect INQUAir Compressor to the nebulizer. Switch it on (For the How-to-use of nebulizer kit, see the IFU of Nebulizer).

10. Cautions

- Read this instruction carefully and thoroughly before use.
- This product should not work continuously for more than 4 hours. To prevent overheating, please turn off the INQUAir Compressor for at least 20 minutes between use.
- The filter cotton needs to be kept dry and clean. The filter cotton should be replaced with a new piece when its color starts to change.
- This product is for aerosol therapy only and is not intended for other purposes.
- Do not use this product with anesthetics.
- Do not pour medicine drop onto the surface of this product, especially near power switch to avoid the risk of short circuit.
- Turn off this product after use.
- Do not put this product onto the surface of hot objects. Keep it where air flow condition is good.
- Do not use this product in humid environment.
- Do not put this product into water or other solutions.
- Do not use this product before sleep to avoid using it unconsciously.
- Children are not allowed to use this product without adult supervision.

Compressor model ANBB26024 is expected to be used under the following electromagnetic environment.

Interference Resistance Test	IEC 60601 Test Electrical Level	Consistent Electrical Level	Electromagnetic Environment - Guideline
Radio Frequency Conductivity on GB/T 17626.6	3 V (Valid Value) 150 KHz ~ 80 MHz	3V (Valid Value)	Portable and mobile radio frequency communication device should not be used closer than the recommended safe distance to any part of the compressor model ANBB26024 including cable. The safe distance is generated by the following formula involving transmitter frequency. Recommended distance d = 80 MHz ~ 800 MHz d = 800 MHz ~ 2.5 GHz —Max. rated power output in watt provided by the transmitter's manufacturer d—Recommended safe distance in meters". Fixed radio frequency transmitter's field intensity is determined by the detection of electromagnetic environment". In each frequency range "should be lower than the consistent electrical level. There may be interference if device with the sign below is nearby.
Radio Frequency Radiation GB/T 17626.3	3V/m 80 MHz ~ 2.5 GHz	3V/m	The power frequency magnetic field should be of the quality that is acceptable by typical commercial or hospital environment

Note 1: at 80MHz and 800MHz frequency, higher frequency range formula is adopted.
Note 2: the guidelines may not be applicable for certain situations. Electromagnetic transmission is affected by the absorption and reflection of buildings, objects and human body.
* The field intensity of fixed transmitter such as wireless phones, ground mobile radio base station, amateur radio, amplitude modulation, frequency modulation broadcast and TV broadcast, cannot be predicted in theory. The survey of the electromagnetic venue should be considered to evaluate the electromagnetic environment of the fixed radio frequency transmitter. If the compressor model ANBB26024 is located in a place where the field intensity is higher than the applicable radio frequency consistent electrical level, observation is needed to verify the functioning of the compressor model ANBB26024. If malfunction is detected, additional measurements must be taken such as adjusting compressor's location or direction.
* In the range between 150KHz and 80MHz, the field intensity should be lower than 3V/m.

10 of 14

12 of 14

13 of 14

14 of 14