
PARI Prüfgerät PG 041

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Stand der Information: September 2011


en Instructions for use

Read these instructions carefully and completely before using for the first time. Keep these instructions in a safe place. Personal safety cannot be guaranteed and/or damage to the device cannot be ruled out if these instructions are not strictly followed.

1 Components

Please check that all components shown are included in your package. If they are not, notify PARI GmbH immediately.

1	Test manometer (0 bar – 2.5 bar; Cl.: 2.5)
2	Test nozzle (Ø 0.48 mm)
3	Connection tubing (1.2 m)
4	Tubing adapter

2 Preparation for use

To prepare the PARI PG 041 test device for use, proceed as follows:

- Insert the reference nozzle in the right side of the T-fitting.
- Insert the connection tubing in the left side of the T-fitting.
- Connect the test device to the compressor using the connection tubing and the tubing adapter.

 If you are working with a PARI COMPACT, use the connection tubing without the adapter.

- Switch the compressor on and read the operating pressure shown on the pressure gauge.

3 Information

The PARI PG 041 test device is an analogue device for measuring the operating pressure of compressors.

The test device is equipped with a precision test nozzle that has a diameter of 0.48 mm. This is the same diameter as that of PARI nebuliser nozzles.

The tolerances listed in the following indicate the pressure range for each compressor type within which PARI nebulisers function optimally. Therefore, it is important and essential to check the compressor pressure regularly to ensure maximum effectiveness of the inhalation treatment.

Important: The PARI PG 041 is a test device. It is therefore recommended to have the device recalibrated three years after it is purchased as new and every two years after that. Please contact PARI GmbH or another authorised Service Center for this.

4 Technical data

Tolerance of the test nozzle at 1000 mbar: +3% / -1%

The minimum pressures that must be reached by the individual PARI compressors are listed in the table:

PARI compressor type	Minimum pressure
PARI BOY SX (Type 085) PARI JuniorBOY SX (Type 085) PARI TurboBOY SX (Type 085)	min. 1.2 bar
PARI BOY mobile S PARI BASIC PARI mini	min. 0.7 bar
PARI SINUS	min. 1.1 bar
PARI COMPACT	min. 0.95 bar

5 Cleaning

The connection tubing is not heat-resistant. It must not be subjected to thermal disinfection or sterilisation procedures. If the connection tubing becomes dirty, it should be replaced as necessary.

The pressure gauge, T-fitting and test nozzle can be wiped off with a damp cloth. Take care to ensure that water does not get into the T-fitting or the test nozzle!

6 Terms and conditions of warranty

During the warranty period (6 months), PARI or our customer service will repair any defects in the PARI PG 041 resulting from faults in material or workmanship free of charge. The owner shall not be entitled to cancel the sale or demand a partial or full refund of the purchase price. This warranty does not cover damage caused by improper use of the PARI PG 041. The warranty shall become void if repairs are carried out on the PARI PG 041 by unauthorised persons. No compensation will be paid for consequential or direct damage.

In the event of a claim, please bring the complete device to your dealer together with your receipt, or send it to PARI Service, packed in the original carton and postage paid.

Wearing parts are excluded from the warranty.

The warranty period shall commence on the date of purchase.

7 Spare parts and accessories

Description	Art. no.
Test nozzle Ø 0.48 mm with O-ring	041B1023
Pressure gauge with T-fitting	041B1030
Connection tubing (1.2 m, f/f)	041B4592
Adapter	041E2802

For all product information please contact our Service Center:
Tel.: +49 (0) 8151-279 279

We reserve the right to make technical changes.

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Pour toute information concernant nos produits, contactez notre S.A.V. : tél. +49 (0)8151-279 279