
# LL Vernebler

CE 0123

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# Hygienic re-use in the professional environment

Information as of: 2012-11. We reserve the right to make technical changes.

<b>GENERAL</b> .....	18
<b>MATERIAL INFORMATION</b> .....	20
<b>PROCESSING LIMITS</b> .....	20
<b>PREPARATION</b> .....	21
<b>CARE OF THE CONNECTION TUBING</b> .....	23
<b>CLEANING AND DISINFECTION</b> .....	24
Mechanical cleaning and disinfection.....	24
Manual cleaning .....	25
Manual disinfection .....	26
<b>STERILISATION</b> .....	27
<b>CLEANING THE NOZZLE</b> .....	29
<b>VISUAL INSPECTION</b> .....	30
<b>STORAGE</b> .....	30
<b>CONTACT</b> .....	30

# 1 GENERAL

---

## NOTICE

Ensure that the system is processed for re-use only with the application of adequately validated device- and product-specific procedures, and that the validated parameters are maintained during each cycle. The equipment used must be maintained regularly.

---

Follow the schedule defined below for hygienic re-use:

– Single patient use:

Clean the nebuliser and the accessories used **after every inhalation session**.

Disinfect the nebuliser and accessories used at least once a day.

– Multiple patient use:

Clean, disinfect and sterilise the nebuliser and the accessories used **before every change in patient**.

Replace the connection tubing because this cannot be processed for re-use by multiple patients.

PARI nebulisers must be processed for re-use with the application of suitable, validated methods in such manner that the success of these methods can be monitored reliably and the safety and health of patients is not threatened.

---

## CAUTION

Ensure that the components are allowed to dry sufficiently after each processing step. Condensation or residual wetness can present an increased risk through the growth of bacteria.

---

The following text describes the methods validated by PARI, with which plastic PARI parts can be processed effectively.

Ensure that the desired results are achieved consistently by the processing methods employed by your staff with the equipment and chemicals used. Take particular care to ensure that the processing procedure you select is effective if it must differ from our validated procedures (e.g., because the equipment and chemicals described are not available). When selecting the chemicals for use, make sure that they are suitable for cleaning or disinfecting medical products made from the materials indicated. Bear in mind the manufacturer's dosing information and instructions for safe handling of chemicals.

To assist you in processing the system for hygienic re-use, a checklist may be obtained from PARI GmbH (on the internet at [www.pari.de/en/products](http://www.pari.de/en/products) or from the PARI Service Center [see also: 11 Contact]).

## 2 MATERIAL INFORMATION

Polypropylene	Inspiratory valve, nebuliser upper section, nebuliser lower section, chamber, medication cup, mouthpiece
Silicone	Valve shim
Polyamide	LL adapter
Fluoroelastomer	Nozzle gasket
Stainless steel	Airflow control, compression spring (for the interrupter button)

## 3 PROCESSING LIMITS

The nebuliser can be sterilised up to 300 times and used for up to two years.


Do not expose the nebuliser or accessories used to a temperature above 137 °C.

The connection tubing cannot be processed for re-use.


## 4 PREPARATION

Dismantle the nebuliser into its individual parts:


- Pull the connection tubing together with the LL adapter off the nebuliser.
- Disconnect the LL adapter from the connection tubing.


- Pull the mouthpiece or mask off the nebuliser.


- Carefully pull the expiratory valve out of the slot in the mouthpiece. The valve must still be attached to the mouthpiece.

- Detach the inspiratory valve and the chamber from the nebuliser.
- Remove the valve shim from the inspiratory valve.

- Unscrew the nebuliser upper section from the nebuliser lower section.
- Push the medication cup up as far as the windows and then pull it out of the nebuliser lower section.
- Press the nozzle together with the airflow control down and out of the medication cup.
- Pull the airflow control up and off from the nozzle.
- Turn the interrupter button clockwise as far as it will go.

## **5 CARE OF THE CONNECTION TUBING**

- If there is condensation visible in the tubing, disconnect the tubing from the nebuliser after the treatment.
- Allow the compressor or the PARI CENTRAL to continue blowing air through the tubing until all of the condensate has been blown out.
- Replace the connection tubing in the case of a change of patients or if it becomes contaminated.

## 6 CLEANING AND DISINFECTION

Clean and disinfect immediately after use. Ideally, a mechanical procedure (instrument dishwasher) is to be used.

---

### NOTICE

When selecting the cleaning agent and/or disinfectant, ensure that it is compatible with the materials used, **otherwise damage to the product components** cannot be ruled out.

---

#### Information when using masks with rubber bands:

**i** The rubber band can only be disinfected chemically, not thermally as it might be damaged by the high temperatures used.

### 6.1 Mechanical cleaning and disinfection

#### EQUIPMENT

- Cleaning and disinfection device in conformance with DIN EN ISO 15883
- Cleaning agent
- Neutralising agent if necessary

#### PROCEDURE

- Arrange the components so that they can be cleaned most effectively.
- Select a suitable cleaning programme and agent.

#### VALIDATION

The method has been validated using:

- G7836 CD cleaner-disinfector produced by Miele, programme Vario TD.
- Standard baskets manufactured by Miele
- Neodisher Mediclean forte 0.5% and Neodisher Z for neutralisation (from Dr. Weigert)

## DRYING

Even if the cleaner-disinfector is equipped with a drying function, make sure that no residual moisture is left in the components. If necessary, remove any residual moisture by shaking the parts and allow them to dry completely.

## 6.2 Manual cleaning

As an alternative to mechanical cleaning and disinfection, the following manual method can be used:

### EQUIPMENT

- Cleaning agent
- Brush
- Warm water

### PROCEDURE

- Place all the parts in the cleaning solution.
- Clean the parts thoroughly with a brush.
- Thoroughly rinse the parts under hygienically tested water and leave them to dry completely on a clean, dry and absorbent surface.

### VALIDATION

The method has been validated using:

- Korsolex® Endo Cleaner 0.5% (Bode)
- Approx. 50 °C warm water
- Duration of cleaning procedure: approx. 5 minutes

## 6.3 Manual disinfection

For manual disinfection, select a suitable instrument disinfectant that has been proven to be effective and compatible with the cleaning agent used.

### EQUIPMENT

- Disinfectant
- Water

---

### CAUTION

PARI product components are not suitable for treatment with instrument disinfectants based on **quaternary ammonium compounds**. Quaternary ammonium compounds tend to become concentrated in plastic materials, which can cause **incompatibility reactions** in patients.

---

### PROCEDURE

- Place all the parts in the disinfection solution.
- Allow the solution to take effect as indicated in the instructions for use of the disinfectant.
- Thoroughly rinse the parts under hygienically tested water and leave them to dry completely on a clean, dry and absorbent surface.

### VALIDATION

The method has been validated using:

- Korsolex® Basic 4% (Bode)
- Treatment time: 15 minutes

## 7 STERILISATION

---


### NOTICE

A validated sterilisation process can only be carried out on products that have been cleaned and disinfected.

---

#### **Information on sterilising the child/adult mask soft:**

Always use the mask stabiliser when sterilising these masks, because otherwise they may lose their shape under the effects of high temperatures. For this purpose, follow the instructions for use of the mask stabiliser as well.

### EQUIPMENT

Steam steriliser – preferably with fractionated pre-vacuum – in accordance with DIN EN 285 or DIN EN 13060 (type B)

### PROCEDURE

- Pack all of the disassembled parts in a sterile barrier system in accordance with DIN EN 11607 (e.g., foil-paper packaging).

Maximum sterilisation temperature and holding time:

121 °C, at least 20 min. or

132 °C / 134 °C, at least 3 min.

## VALIDATION



The method has been validated according to  
DIN EN ISO 17665-1 using:

- Euro-Selectomat 666 from MMM (fractionated pre-vacuum method)
- Varioklav 400 E from H+P Labortechnik (gravitation method)

## 8 CLEANING THE NOZZLE

Should the nozzle be blocked, you can clear it yourself as follows:

- Unscrew the nebuliser upper section from the nebuliser lower section.
- Push the medication cup up as far as the windows and then pull it out of the nebuliser lower section.
- Press the nozzle (5) together with the airflow control (4) down and out of the medication cup.
- Pull the airflow control up and off from the nozzle.
- Take the nozzle cleaner (6) out of the nebuliser lower section (1).
- Carefully poke through the three nozzle bore holes "A" to "C" (only from below).
- Reassemble the nebuliser.

### NOTICE

There is a **risk of damaging the nozzle bore holes irreparably**. Therefore, always use only an original PARI nozzle cleaner and proceed extremely carefully. If the nozzle becomes blocked repeatedly, replace the nozzle.

## **9 VISUAL INSPECTION**

Inspect all product components after each cleaning, disinfection and sterilisation. Replace any broken, misshapen or seriously discoloured parts.

## **10 STORAGE**

Store the nebuliser and the accessories used in a dry, dust-free location equipped with protection against contamination. It is recommended to use sterile packaging.

## **11 CONTACT**

For all product information and in the event of defects or questions about usage, please contact our Service Center:

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