
## **VORTEX® TRACHEO**

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Information as of: April 2017. All rights reserved. Technical and design changes and printing errors reserved. Principle shown in illustrations.

<b>IMPORTANT INFORMATION</b> .....	<b>29</b>
General.....	29
Information about the instructions for use .....	29
Structure of safety instructions .....	29
Treatment of babies, children and anyone who requires assistance .....	30
Hygiene .....	30
<b>PRODUCT DESCRIPTION</b> .....	<b>31</b>
Components .....	31
Intended purpose .....	31
Intended use.....	32
Contraindications.....	32
Description of function .....	32
Material information .....	33
Operating life .....	33
<b>INHALATION</b> .....	<b>34</b>
Function check .....	34
Preparing for inhalation .....	35
Performing the inhalation .....	36
<b>CLEANING AND DISINFECTION AT HOME</b> .....	<b>39</b>
Preparation .....	39
Cleaning .....	40
Disinfection .....	40
Visual inspection .....	41
Drying and storage .....	41

<b>HYGIENIC RE-USE IN THE PROFESSIONAL ENVIRONMENT .....</b>	<b>42</b>
General.....	42
Processing limits .....	43
Preparation .....	43
Cleaning and disinfection .....	44
Sterilisation .....	47
Visual inspection .....	48
Storage .....	48
<b>MISCELLANEOUS .....</b>	<b>49</b>
Disposal.....	49
Explanation of symbols .....	49
Contact .....	49

# 1 IMPORTANT INFORMATION

## 1.1 General

Please read these instructions for use carefully and completely. Do not discard, so you can consult them at a later date. If you fail to comply with the instructions for use, injury or damage to the device cannot be ruled out.

**If symptoms persist or if your condition worsens, discontinue the application and contact your doctor immediately.**

## 1.2 Information about the instructions for use

If these instructions for use are lost, you can request another copy from PARI GmbH [see: Contact, page 49]. The instructions for use for some products can be retrieved on the internet in both German and English. Simply visit: [www.pari.com](http://www.pari.com) (on the respective product page).

## 1.3 Structure of safety instructions

Safety-critical warnings are categorised according to hazard levels in these instructions for use:

- The signal word **WARNING** is used to indicate hazards which, without precautionary measures, can result in serious injury or even death.
- The signal word **CAUTION** is used to indicate hazards which, without precautionary measures, can result in minor to moderate injury or impair treatment.
- The signal word **NOTICE** is used to indicate general precautionary measures which are to be observed to avoid damaging the product during use.

## **1.4 Treatment of babies, children and anyone who requires assistance**

Babies, children and anyone who requires assistance must be supervised constantly by an adult during inhalation therapy. This is the only way to ensure safe and effective treatment. Individuals in this group often underestimate the hazards involved, thus resulting in a risk of injury.

The product includes small parts, which can pose a choking hazard. Therefore, keep this product out of the reach of children younger than 3 years.

This product is only designed for patients who are able to breathe by themselves and are conscious. This is the only way to ensure effective treatment and avoid a choking hazard.

## **1.5 Hygiene**


Please observe the following hygiene instructions:

- When used in the home environment:  
Each VORTEX must only be used by a single patient for hygiene reasons.
- When used in the medical environment:  
Please observe the general hygiene regulations that apply in doctors' practices/hospitals.
- Make absolutely sure you also carry out cleaning and drying or hygienic preparation for re-use before using for the first time.
- For cleaning and disinfection, always use fresh tap water of drinking water quality.
- Make sure all components are allowed to dry properly after every cleaning or disinfection.
- Do not keep the VORTEX and accessories in a damp environment or together with damp objects.

## 2 PRODUCT DESCRIPTION

### 2.1 Components

Check that all components of your PARI product are contained in your package. If anything is missing, please notify the dealer from whom you purchased the PARI product immediately.

- (1) Inhalation chamber (antistatic)
  - (1a) Connecting ring for metered dose inhaler (not removable)
- (2) Mouthpiece with inspiratory and expiratory valve
  - (2a) Protective cap
- (3) Tracheal adapter

### 2.2 Intended purpose

The VORTEX is a holding chamber, and is used in conjunction with medication sprays or "metered dose inhalers" in the treatment of diseases of the airways.

The VORTEX Tracheo is equipped with a tracheal adapter. In conjunction with the tracheal adapter, the VORTEX is suitable for use by tracheostomised patients who are not able to inhale through a mouthpiece.

The tracheal adapter is suitable for use with blocked (with cuff) and unblocked tracheal cannulas. However, it must not be used on an open tracheostomy.

## **2.3 Intended use**

The VORTEX is a holding chamber which makes it easier to inhale from a metered dose inhaler.

The VORTEX Tracheo must not be used by tracheostomised patients without professional instruction. It can also be used for babies and children under constant adult supervision.

Only medication that has been approved for inhalation treatment must be used. Take note of any restrictions in the instructions for use of the medication in question.

In order to ensure proper treatment, please use only PARI accessories.

## **2.4 Contraindications**

There are no contraindications known to PARI GmbH.

## **2.5 Description of function**

The VORTEX with its metal chamber and cyclone twist principle guarantees that medication is dispensed reliably in metered quantities. It is used to minimise coordination mistakes when using metered dose inhalers.

With a flexible connecting ring it can be used with all standard metered dose inhalers.

## 2.6 Material information

Aluminium	Inhalation chamber
Polymethylpentene	Mouthpiece
Polypropylene	Connecting ring for metered dose inhaler, mouthpiece, tracheal adapter, protective cap
Silicone	Inspiratory and expiratory valve in the mouthpiece
Thermoplastic elastomer	Connecting ring for metered dose inhaler, mouthpiece

All materials are latex-free.

## 2.7 Operating life

Replace the VORTEX regularly according to usage:

### **When used in the home environment**

The VORTEX should be replaced after 60 disinfection cycles, if used daily after not more than a year.

### **When used in the medical environment**

The VORTEX should be replaced after it has been prepared hygienically for re-use 30 times, at the latest after not more than a year.


## 3 INHALATION

### 3.1 Function check

#### Clear pathway through the valve

**After receiving the VORTEX**, check that the path through the valve in the mouthpiece is unobstructed:


- Run tap water through the lip attachment into the mouthpiece.  
The water should escape through the side openings in the mouthpiece.
  - Let the tap water run out through the bottom of the mouthpiece.  
The water should run out of the lip attachment again.
  - Dry the VORTEX completely before using it for the first time.
- If the valve is obstructed, the VORTEX must not be used.

#### Position of the valve

**Before every use**, check that the valve is in the correct position inside the VORTEX mouthpiece:

- Check that the valve is in the position shown inside the mouthpiece.  
If necessary, use a stream of water to correct its position.
- Dry the VORTEX completely before using it again.

## 3.2 Preparing for inhalation

### WARNING

Since the inhalation chamber of the VORTEX is not completely closed, small particles can get in, and might be breathed in during an inhalation session (**choking hazard**). Therefore, it is imperative to make sure there are no foreign bodies inside the VORTEX before every use.

### CAUTION

Check all product components and the accessories before each use. Replace any broken, misshapen or seriously discoloured parts. Please also follow the instructions for assembly given below. Damaged components and an incorrectly assembled holding chamber may impair the function of the holding chamber and thus the **treatment as well**.


### CAUTION

**Risk of infection** in the case of a change in patients. Always use a new tracheal adapter when changing patients. If the VORTEX was not sterilised before a change in patients, a new mouthpiece must also be used.


- Attach the mouthpiece to the inhalation chamber.

- Remove the protective cap on the metered dose inhaler.
- Shake the metered dose inhaler vigorously before every use.

- Insert the metered dose inhaler in the connecting ring on the inhalation chamber.
- Take the protective cap off the mouthpiece.

- Place the tracheal adapter on the mouthpiece.

### 3.3 Performing the inhalation

#### CAUTION


Before the patient starts inhaling, ensure that all parts are firmly connected to each other. Otherwise, **insufficient** dosage may be delivered.

**⚠ CAUTION**


Children younger than 18 months should inhale four to six times; children older than 18 months should inhale deeply two to four times. Otherwise, insufficient **dosage may be delivered**.

For adults, one breath is enough to inhale the medication from the VORTEX.


***A useful tip for visually checking breathing (e.g., with infants): The valves must visibly open when exhaling.***


- Assume a posture in which the upper body is upright (sitting or lying with the upper body raised).
- Connect the tracheal adapter and the tracheal cannula.


- Hold the VORTEX as shown in the illustration.

- Actuate the metered dose inhaler at the start of a deep, slow intake of breath.
- Hold your breath briefly.


- Breathe out slowly into the VORTEX.

- Disconnect the tracheal adapter together with the VORTEX from the tracheal cannula.

**⚠ CAUTION! Risk of injuring the tracheostomy!** Hold the tracheal cannula firmly as you pull the tracheal adapter out.

## 4 CLEANING AND DISINFECTION AT HOME

Regularly clean all product components and the accessories used, if they are visibly dirty immediately after use.

Disinfect all individual parts at least once a week, in the event of an acute infection, once a day.

### 4.1 Preparation

- Take the metered dose inhaler out of the connecting ring on the inhalation chamber and close it with the protective cap provided for this purpose.
- Remove the tracheal adapter from the mouthpiece.
- Disconnect the mouthpiece from the inhalation chamber.

**Info:** *The protective cap must be open.*

## 4.2 Cleaning

- Place all disassembled components in warm tap water with a little dishwashing liquid for at least 5 min.
- Rinse all parts thoroughly in running water.
- You can remove excess water more quickly by shaking all parts.

## 4.3 Disinfection

Once they have been cleaned, disinfect the VORTEX and accessories as described in the following.

### CAUTION

A damp environment may encourage the growth of bacteria. Therefore, remove all parts from the pot or disinfectant as soon as disinfection has finished. Dry the parts. The **risk of infection** is reduced when the parts are dried completely.

## In boiling water

- Place all the **individual parts** in boiling water for at least 5 minutes. Use a clean pot and fresh water with reduced calcium content, if possible.

**NOTICE! Risk of damaging plastic parts!** Plastic will melt if it comes into contact with the hot base of the pot. Make sure the water level in the pot is high enough to enable the parts to float clear of the base.

- You can remove excess water more quickly by shaking all parts.

## Using a standard thermal disinfecter for baby bottles (not a microwave oven)

For effective disinfection, use a disinfecter with a runtime of at least 6 minutes. Regarding disinfection, the duration of the disinfection procedure and the quantity of water required for this, follow the instructions for use of the disinfecter you are using.

### CAUTION

Inadequate disinfection encourages the growth of bacteria and thus increases the **risk of infection**. Thorough disinfection has not been completed until the disinfecter automatically switches itself off, or the minimum disinfection time specified in the instructions for use of the disinfecter has elapsed. Therefore, do not switch the device off prematurely. Also make sure that the disinfecter is kept clean and regularly check that it is in good working order.

## 4.4 Visual inspection

Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discoloured parts.

## 4.5 Drying and storage

- Dry all individual parts completely.  
Because of the antistatic properties of the inhalation chamber, you can dry them with a clean tea towel.
- Place the protective cap over the mouthpiece.  
**Info:** *The protective cap protects the mouthpiece from damage during transport and storage.*
- Attach the mouthpiece to the inhalation chamber.
- Store the VORTEX and its accessories in a dry, dust-free place (safely out of continuous direct sunlight).

## 5 HYGIENIC RE-USE IN THE PROFESSIONAL ENVIRONMENT

### 5.1 General

#### CAUTION

**Risk of infection** in the case of a change in patients. The VORTEX must be cleaned, disinfected **and** sterilised before a change of patients. The tracheal adapter must be disposed of. If sterilisation is not carried out, the mouthpiece must also be disposed of.

#### CAUTION

Ensure that the components are allowed to dry sufficiently after each processing step. Condensation or residual wetness can present an increased **risk of bacterial growth**.

#### **NOTICE**

Ensure that the system is processed for re-use only with the application of adequately validated device- and product-specific procedures, and that the validated parameters are maintained during each cycle. The equipment used must be maintained regularly.

PARI products must be processed for re-use with the application of suitable, validated methods in such manner that the success of these methods can be monitored reliably and the safety and health of patients is not threatened.

The following text describes the methods validated by PARI, with which the product components can be processed effectively. If you have any questions about processing for re-use, please contact the PARI Service Center [see: Contact, page 49].

## Processing cycles

Observe the following cycle for hygienic processing of the VORTEX and the accessories used:

– **Single patient use:**

Clean and disinfect all product components and the accessories used

- once a day in the case of acute infection of the airways,
- at least once a week,
- whenever there is visible soiling.

– **Multiple patient use:**

Clean, disinfect and sterilise all product components before every change in patient.

Sterilisation may be omitted if a new mouthpiece is used.

## 5.2 Processing limits

The VORTEX can be processed for re-use up to 30 times and used for up to one year.

## 5.3 Preparation

Disassemble the VORTEX:

- Take the metered dose inhaler out of the connecting ring on the inhalation chamber and close it with the protective cap provided for this purpose.
- Remove the tracheal adapter from the mouthpiece.
- Disconnect the mouthpiece from the inhalation chamber.
- In the event of a change of patients without prior sterilisation of the VORTEX, dispose of the mouthpiece and the tracheal adapter.

**Info:** *The protective cap must be open.*

- First rinse all components briefly under running tap water. This must be done especially thoroughly if cleaning and disinfection are not carried out straight away.

## 5.4 Cleaning and disinfection

Generally, a mechanical process (instrument dishwasher) should be used for cleaning and disinfection.

### NOTICE

When selecting the cleaning agent and/or disinfectant, ensure that it is compatible with the materials used, **otherwise damage to the product components** cannot be ruled out.

## Mechanical cleaning and disinfection

### EQUIPMENT

- Cleaning and disinfection device in conformance with DIN EN ISO 15883
- Cleaning agent
- Neutralising agent if necessary

### PROCEDURE

- Arrange the **components** so that they can be cleaned most effectively.  
The inhalation chamber must be held upright.
- Select a suitable cleaning programme and agent.

**i** *In rare cases, cleaning in the instrument dishwasher can cause discolouration of the aluminium. This is harmless and does not impair the functioning of the VORTEX.*

## VALIDATION

The method has been validated in Europe using:

- RDG G7836 CD cleaner-disinfector produced by Miele (programme Vario TD)<sup>1</sup>
- Miele standard baskets (E450 / E142)
- Neodisher MediZym as pH-neutral cleaning agent (0.5%, by Dr. Weigert)

**Alternatively:** Neodisher Mediclean forte 0.5% (alkaline), and for neutralisation Neodisher Z (by Dr. Weigert)

## DRYING

Even if the cleaner-disinfector is equipped with a drying function, make sure that no residual moisture is left in the components. If necessary, remove any residual moisture by shaking the parts and allow them to dry completely.

## Chemical manual cleaning and disinfection

Alternatively, cleaning and disinfection can also be carried out manually. The validated agents are listed together with all relevant information in the following table:

	Korsolex® Endo Cleaner 0.5% (Bode)	Korsolex® Basic 4% (Bode)
Cleaning agent	×	–
Disinfectant	–	×
Contains aldehyde?	formaldehyde-free	×
Concentration	0.5%	4%
Duration	5 min.	30 min.
Temperature	50 °C	max. room temperature

<sup>1</sup>) Cleaning phase at temperatures up to 55 °C: 5 minutes' holding time.  
Thermal disinfection at 90 °C with 5 minutes' holding time.

## EQUIPMENT

- Cleaning agent/Disinfectant
- Warm water
- Brush if applicable
- Basin

## PROCEDURE

### CAUTION

Inadequate cleaning and disinfection encourages the growth of bacteria and thus increases the **risk of infection**. Adequate cleaning with disinfection can only be assured if the specified mixing ratio and application time are adhered to, and if all **individual parts** are completely immersed in the solution for the entire application time. There must not be any air pockets or bubbles.

### CAUTION

PARI product components are not suitable for treatment with instrument disinfectants based on **quaternary ammonium compounds**. Quaternary ammonium compounds tend to become concentrated in plastic materials, which can cause **incompatibility reactions** in patients.

If a disinfectant without a cleaning function is used, the following procedure must be performed twice: the first time with a cleaning agent, and then with the disinfectant.

To ensure safety when handling chemicals, follow the instructions for use of the agent used, particularly the accompanying safety instructions.

- Prepare the solution in the basin in accordance with the manufacturer's instructions. The quantity of solution required depends on the number of parts that are to be cleaned.

- Place all the individual parts in the prepared solution and leave them to soak.  
*Info: If the application period is exceeded significantly, the plastic parts may take on the smell of the disinfectant.*
- Rinse off all parts thoroughly in running water (residues of the disinfectant can cause allergic reactions or irritations of the mucous membrane).
- Leave all parts to dry completely on a dry, clean, absorbent surface.
- Dispose of the used solution.  
The diluted solution can be disposed of down the drain.

## 5.5 Sterilisation

### NOTICE

A validated sterilisation process can only be carried out on products that have been cleaned and disinfected.

 *The tracheal adapter is not sterilisable.*

### EQUIPMENT

Steam steriliser – preferably with fractionated pre-vacuum – in accordance with DIN EN 285 or DIN EN 13060 (type B)

### PROCEDURE

- Pack all of the disassembled parts in a sterile barrier system in accordance with DIN EN 11607 (e.g., foil-paper packaging).

Maximum sterilisation temperature and holding time:

121 °C, at least 20 min. or

132 °C / 134 °C, at least 3 min.

## VALIDATION

The method has been validated according to DIN EN ISO 17665-1 using:

- Euro-Selectomat 666 from MMM (fractionated pre-vacuum method)
- Varioklav 400 E from H+P Labortechnik (gravitation method)

## 5.6 Visual inspection

Inspect all product components after each cleaning, disinfection and sterilisation. Replace any broken, misshapen or seriously discoloured parts.

## 5.7 Storage

Store the assembled VORTEX and the accessories used in a dry, dust-free location equipped with protection against contamination.

## 6 MISCELLANEOUS

### 6.1 Disposal

All product components can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.

### 6.2 Explanation of symbols

The following symbols are displayed on the packaging of your PARI product:

	Manufacturer
	Please follow the instructions for use.
	The product satisfies the basic requirements as set forth in Appendix I of Directive 93/42/EEC concerning medical devices.
	Order no.
	Lot no.

### 6.3 Contact

For all product information and in the event of defects or questions about usage, please contact our Service Center:

Tel.: +49 (0)8151-279 279 (German-speaking)  
+49 (0)8151-279 220 (international)

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