


Zirela® nebuliser handset

REF 678G8015

de	Gebrauchsanweisung	fr	Notice d'utilisation
en	Instructions for Use	it	Istruzioni per l'uso
da	Brugsanvisning	nl	Gebruiksaanwijzing
es	Instrucciones de uso	sv	Bruksanvisning
fi	Käyttöohje		

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Zirela® Vernebler

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en Instructions for use

Zirela® nebuliser handset

- *Read these instructions for use carefully and completely.*
- *Read the **instructions for use of the eFlow®rapid nebuliser system completely before using for the first time. These contain further important instructions for the use of the Zirela® nebuliser handset, particularly regarding connection, operation, cleaning and disinfection.***

1 Intended purpose

The Zirela® nebuliser handset with an eBase® Controller is intended to administer Quinsair® as an aerosol for inhalation.

2 Indication

The Zirela® nebuliser handset is indicated to administer Quinsair® to treat lung infections caused by *Pseudomonas aeruginosa* in patients with cystic fibrosis (CF).

3 Contraindication

For the Zirela® nebuliser handset no contraindication is known to PARI Pharma GmbH. Take note of contraindications of the medicine in the Quinsair® package leaflet.

4 Use of the nebuliser handset

The Zirela® nebuliser handset must only be used in accordance with its intended purpose. For reasons of hygiene and to avoid cross-infections, the Zirela® nebuliser handset must only be used by one patient.

- Before using for the first time, disassemble and clean the nebuliser handset including the aerosol head as described in the instructions for use of the eFlow®rapid nebuliser system.
- Assemble the nebuliser handset as shown in the diagram:
 - Open the nebuliser handset.
 - Insert the aerosol head. **Do not touch the membrane (the curved metal surface in the middle of the aerosol head).**
 - Insert the inspiratory valve.
 - Close the nebuliser handset.
 - Attach the mouthpiece.
 - Connect the connection cord.
 - Fill with the medication; empty an ampoule completely.
 - Close the medication reservoir.
 - Sit in an upright position and relax. This makes inhaling easier and improves uptake of the medication in the airways.

-  Hold the nebuliser handset horizontally all the time while performing the inhalation. If the handset is tilted, errors in the dosage may occur.

- Hold the mouthpiece between your teeth and enclose it with your lips. The lips must not touch the blue expiratory valve.

The inhalation session should last about 5 minutes. Inhalation is ended as soon as no more aerosol is generated.

- Always check the medication reservoir to ensure that the entire dose of Quinsair® has been nebulised.

5 Cleaning and disinfection

The dismantled nebuliser handset (including the aerosol head) must be cleaned immediately after each use and disinfected at least once a day (additional information can be found in the instructions for use of the eFlow®rapid nebuliser system).

6 Replacing the nebuliser handset

Each 28-day pack of Quinsair® contains a Zirela® nebuliser handset. Use the new nebuliser handset for each treatment cycle.

7 Package contents

One Zirela® nebuliser handset with one Zirela® aerosol head. There are no accessories for the Zirela® nebuliser handset.

8 Symbols

	Medical device
	Consult instructions for use
	Single patient – Multiple use
	Batch code
	Item number
	Serial number
	Minimum and maximum ambient temperature
	Minimum and maximum atmospheric humidity
	Minimum and maximum atmospheric pressure
	CE marking: This product is compliant with the requirements set forth in Directive 93/42/EEC (Medical Devices).
	The medical device was distributed commercially after 13 August 2005. This product must not be disposed of with normal domestic waste. The symbol of a crossed-through refuse bin indicates that separate collection is required.
	Caution
	Manufacturer
	Date of manufacture

9 Service and contact

In case of technical problems or questions about the device, please contact your service partner. The contact addresses for your responsible service partners as well as a digital version of your Zirela® instructions for use can be found at **www.pari.com/zirela**.

Aerosol data is available from PARI Pharma upon request.

We reserve the right to make technical changes.

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