de	Gebrauchsanweisung	3
en	Instructions for use	37

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en Instructions for use

Information as of: 2015-05. All rights reserved. Technical and design changes and printing errors reserved. Principle shown in illustrations.

IMPORTANT INFORMATION	39
General	39
Information about the instructions for use	39 39
Treatment of babies, children and anyone who requires assistance	40
Hygiene	40
PRODUCT DESCRIPTION	41
Components	41
Intended use	42
Proper use	42
Contraindications	43
Product combinations	43
Material information	43 43
Operating life	43
INHALATION	44
Assembling the nebuliser	44
Filling with medication	47
Performing the inhalation	48
CLEANING AND DISINFECTION AT HOME	51
Preparation	51
Care of the connection tubing	51
Cleaning	52
Disinfection	52
Visual inspection	57
Drying and storage	57

HYGIENIC RE-USE IN THE PROFESSIONAL ENVIRONMENT	57
General	57
Processing limits	58
Preparation	59
Cleaning and disinfection	59
Sterilisation	64
Visual inspection	65
Storage	65
MISCELLANEOUS	65
Disposal	65
Technical data	65
Aerosol characteristics	66
Explanation of symbols	68
Contact	68

1 IMPORTANT INFORMATION

1.1 General

Please read these instructions for use carefully and completely. Do not discard, so you can consult them at a later date. If you fail to comply with the instructions for use, injury or damage to the device cannot be ruled out.

If symptoms persist or if your condition worsens, discontinue the application and contact your doctor immediately.

1.2 Information about the instructions for use

If these instructions for use are lost, you can request another copy from PARI GmbH [see: Contact, page 68]. The instructions for use for some products can be retrieved on the internet in both German and English. Simply visit: www.pari.de/de-en/products (on the respective product page).

1.3 Structure of safety instructions

Safety-critical warnings are categorised according to hazard levels in these instructions for use:

- The signal word WARNING is used to indicate hazards which, without precautionary measures, can result in serious injury or even death.
- The signal word CAUTION is used to indicate hazards which, without precautionary measures, can result in minor to moderate injury or impair treatment.
- The signal word NOTICE is used to indicate general precautionary measures which are to be observed to avoid damaging the product during use.

1.4 Treatment of babies, children and anyone who requires assistance

Babies, children and anyone who requires assistance must be supervised constantly by an adult during inhalation therapy. This is the only way to ensure safe and effective treatment. Individuals in this group often underestimate the hazards involved (e.g. strangulation with the power cord or the connection tubing), thus resulting in a risk of injury.

The product contains small parts. Small parts can block the airways and lead to a choking hazard. Therefore, make sure that you always keep the compressor, the nebuliser and the accessories out of the reach of babies and infants.

Special masks may be obtained for treating babies and infants who are not yet able to inhale using the mouthpiece.

The nebulisers are only designed for patients who are able to breathe by themselves and are conscious. This is the only way to ensure effective treatment and avoid a choking hazard.

1.5 Hygiene

Observe the following hygiene instructions:


- Do not use product components unless they have been thoroughly cleaned and dried. Contamination and residual moisture encourage the growth of bacteria, so increasing the risk of infection.
- Before every use and cleaning, wash your hands thoroughly.
- Make absolutely sure you also carry out cleaning and drying before using the device for the first time.
- For cleaning and disinfection, always use fresh water with reduced calcium content.
- Make sure all components are dried properly whenever you have performed cleaning or disinfection.
- Do not keep the product and accessories in a damp environment or together with damp objects.

2 PRODUCT DESCRIPTION

2.1 Components

The components included vary depending on the pack. To find out which components are included in your pack, refer to the accompanying leaflet.

Check that all components of your PARI product are contained in your pack. If anything is missing, please notify the dealer from whom you purchased the PARI product immediately.

- (1) Nebuliser upper section
- (2) Nozzle attachment
- (3) Nebuliser lower section
- (4) Mouthpiece
- (5) Child or adult mask soft
- (6) Baby mask with baby bend
- (7) a) Connection tubing f/m
 - b) Connection tubing f/f with adapter

2.2 Intended use

The PARI XLent nebuliser is a medication nebuliser that is intended for use with a PARI compressor for inhalation therapy of the upper and central airways.

An application usually takes about 5 minutes (depending on the quantity of medication).

Only original PARI compressors and PARI accessories must be used. PARI GmbH cannot be held liable for damage caused by using PARI nebulisers with compressors or accessories from other manufacturers.

Please read and follow the instructions for use of the compressor and of the accessories you are using.

Extended range of intended use for a professional setting

The nebuliser can be operated in conjunction with a PARI compressor or with the PARI CENTRAL¹ in a medical facility. It can be used for multiple patients provided the special requirements regarding hygienic re-use are strictly observed.

2.3 Proper use

For reasons of hygiene, the nebuliser, the mouthpiece and the masks must only be used by a single patient.

This product is only designed for patients who are able to breathe by themselves and are conscious.

Only medication that has been approved for inhalation treatment must be used. Take note of any restrictions in the instructions for use of the medication in question.

The PARI CENTRAL is not available in all countries. Information about the PARI CENTRAL can be obtained from PARI GmbH.

2.4 Contraindications

None.

2.5 Product combinations

The PARI LC SPRINT XLent nebuliser can be operated with all PARI compressors that are available on the market.

To ensure that the nebuliser can be used for all age groups, it is compatible with the corresponding accessories:

- with baby mask and baby bend for treatment of the airways in babies (including premature babies) and infants.
- with mouthpiece or adult mask for general treatment of the airways in children aged approx. 4 years and older and adults.

2.6 Material information

Polypropylene	Nebuliser, nozzle attachment, mouthpiece, baby bend, child mask, adult mask
Silicone	Baby mask
Thermoplastic elastomer	Nebuliser, nozzle attachment, mouthpiece, child mask, adult mask
Polyester	Rubber band
Synthetic rubber (latex-free)	Rubber band

2.7 Operating life

Replace the nebuliser after a year at the latest as changes can occur in the material of the nebuliser over time. This might impair the effectiveness of the treatment.

3 INHALATION

Risk of cross-contamination in the case of a change in patients (in the professional environment): If the nebuliser is used for more than one patient, there is a risk that the germs can be transmitted from one patient to the next via the nebuliser and its accessories. To prevent this, the nebuliser and accessories must be cleaned, disinfected and sterilised, and the connection tubing must be cleaned and disinfected mechanically or replaced before every change of patient [see: General, page 57].

3.1 Assembling the nebuliser

A CAUTION

Check all nebuliser components and the accessories before each use. Replace any broken, misshapen or seriously discoloured parts. Please also follow the instructions for assembly given below. Damaged components and/or an incorrectly assembled nebuliser may impair functioning of the nebuliser and thus treatment as well.

- Press the nozzle attachment lightly onto the nozzle in the nebuliser lower section.
 The arrow on the nozzle attachment must point upwards.
- Place the nebuliser upper section on the nebuliser lower section and turn it clockwise to lock the two parts together.

 Attach the connection tubing to the nebuliser.

Using with a mouthpiece

 Fit the mouthpiece onto the nebuliser.

Use with a baby mask and baby bend

- Attach the baby mask to the baby bend.
 Use your thumb to press the baby mask from the inside over the first ridge on the baby bend.
- Attach the baby bend to the nebuliser.

Using with a child or adult mask

 Attach the child or adult mask to the nebuliser.

When using a mask, make sure that the expiratory valve plate is pressed out so you can exhale freely during an inhalation session.

3.2 Filling with medication

For optimum aerosol generation, fill the nebuliser with at least 3 ml medication.

NOTICE

Take care to ensure that the cap only moves in the direction allowed by the hinge. Otherwise the **nebuliser cap might snap off.**

- Insert the nebuliser in the holder on the compressor intended for this purpose.
- Open the nebuliser cap by pressing your thumb against the underside of the cap.
- Pour the quantity of medication prescribed by your doctor into the nebuliser nozzle attachment.

Be sure to follow the instructions regarding the minimum and maximum fill volumes [see: Technical data, page 65]. If the nebuliser contains too little or too much medication, nebulisation, and thus also the treatment, will be greatly impaired.

Close the nebuliser cap.
 Make sure that the cap snaps into place.

If you are to inhale several medications one after the other in one application, rinse out the nebuliser with tap water between inhalations. Remove excess water by shaking the nebuliser and then fill with medication as described.

3.3 Performing the inhalation

Before the patient starts inhaling, ensure that all parts are firmly connected to each other.

Perform inhalation until the medication is used up (indicated by a change in the sound of the nebuliser).

The nebuliser produces relatively large droplets, and these can be deposited on the nebuliser wall. In order to minimise loss of medication, tap the nebuliser lightly with your finger a few times towards the end of the inhalation session. The droplets will run back down into the nozzle attachment, where they can be nebulised again.

It is normal for a few drops of medication to be left in the nebuliser after the end of the inhalation session.

When using a PARI compressor:

 Connect the connection tubing to the compressor and switch the compressor on.

When using with a PARI CENTRAL:

Prerequisite: A PARI CENTRAL must be connected to the wall outlet fitting in the hospital.

 Insert the connection tubing into the PARI CENTRAL using the adapter provided for this purpose. Twist slightly to make the connection easier.

⚠ WARNING

Check carefully to ensure that the other end of the PARI connection tubing that is connected to the compressor or the PARI CENTRAL is attached to the PARI nebuliser. When several tubing systems are present, there is a possible **threat to life** if different connection options are confused with each other. This applies particularly **in the case of patients in need of constant care** who must receive artificial nutrition or other infusion treatments.

Inhalation with mouthpiece or with child or adult mask

- Sit in an upright position and relax.
- · Take the nebuliser out of its holder and hold it upright.
- Take the mouthpiece between your teeth and close your lips around it or press the mask lightly over the mouth and nose.
- Breathe in through the mouthpiece or mask as slowly and deeply as possible, and then breathe out in a relaxed manner.
- Replace the nebuliser in the holder as soon as you have finished inhaling.
- · Switch the compressor off.

Inhalation with a baby mask and baby bend

- The baby mask must only be used in conjunction with the baby bend. The baby or infant can breathe out without obstruction and the risk of shortness of breath is avoided.
- Adjust the alignment of the bend and the mask according to the position the baby or infant is in.
- Gently press the mask against the face so that it fits snugly over the mouth and nose.
 Make sure that the nebuliser is upright.

- Sick infants often struggle when the mask is pressed against their face, and twist their head back and forth. To ensure effective inhalation, hold the mask firmly in your hand and rest your little finger against the child's cheek. This will enable you to follow the movements of the child's head with the mask more easily.
- · Complete the treatment session.
- Replace the nebuliser in the holder as soon as you have finished inhaling.
- · Switch the compressor off.

4 CLEANING AND DISINFECTION AT HOME

The nebuliser and the accessories used must be cleaned thoroughly after each application and disinfected at least once a day.

The connection tubing cannot be cleaned or disinfected.

4.1 Preparation

- · Detach the tubing from the nebuliser.
- Make sure that all medication residues are removed from the nebuliser.
- Dismantle the nebuliser into its individual parts.
- Carefully pull the blue expiratory valve out of the slot in the mouthpiece. The valve must still be attached to the mouthpiece.


Rinse all components of the nebuliser under running tap water. This must be done especially thoroughly if cleaning and disinfection are not carried out straight away.

4.2 Care of the connection tubing

- · Connect the connection tubing to the compressor.
- Switch the compressor on.
- Leave the compressor running until any condensation in the tubing has been removed by the air flowing into the tube.

4.3 Cleaning

Place all disassembled components in warm tap water with a little dishwashing liquid for at least 5 min. If necessary, use a clean brush to remove loose dirt (the brush must be reserved exclusively for this purpose).

- · Rinse all parts thoroughly in running water.
- You can remove excess water more quickly by shaking all parts.

4.4 Disinfection

After cleaning, disinfect the dismantled nebuliser and accessories used (except for the connection tubing) (effective disinfection is only possible if the nebuliser has been cleaned).

Information when using masks with rubber bands:
The rubber band can only be cleaned as it might be damaged by the high temperatures used during disinfection.

A CAUTION

A damp environment may encourage the growth of bacteria. Therefore, remove all parts from the pot or disinfector as soon as disinfection has finished. Dry the parts. The **risk of infection** is reduced when the parts are dried completely.

In boiling water

- Place all the individual parts in boiling water for at least 5 minutes.
 - Use a clean pot and fresh water with reduced calcium content, if possible.
- Plastic will melt if it comes into contact with the hot base of the pot. Therefore, make sure there is plenty of water in the pot. This way you will avoid damaging the components.
- You can remove excess water more quickly by shaking all parts.

Using a standard thermal disinfector for baby bottles (not a microwave oven)

For effective disinfection, use a disinfector with a runtime of at least 6 minutes. Regarding disinfection, the duration of the disinfection procedure and the quantity of water required for this, follow the instructions for use of the disinfector you are using.

⚠ CAUTION

Inadequate disinfection encourages the growth of bacteria and thus increases the **risk of infection**. Thorough disinfection has not been completed until the disinfector automatically switches itself off, or the minimum disinfection time specified in the instructions for use of the disinfector has elapsed. Therefore, do not switch the device off prematurely. Also make sure that the disinfector is kept clean and regularly check that it is in good working order.

With a microwave disinfector

Microwave disinfectors ensure effective disinfection of the product components in a standard commercial microwave oven.

Place the product components in the disinfector with enough space between them to allow the steam to reach all surfaces reliably.

Maintain the following times when disinfecting:

Wattage	Disinfection time	Cooling time
850 – 950	4 min.	2 min.
600 - 850	6 min.	2 min.

Information about the amount of water required and about the disinfection process is included in the instructions for use of the devices in question (disinfector and microwave).

CAUTION

The product components must not be exposed to microwave radiation unless they are in a microwave disinfector. If they are heated in a microwave oven without the protection of the disinfector, the **individual parts will be damaged.**

Do not use the "grill" function because this too will **damage the individual parts.**

At the end of the disinfection process, the parts may still be very hot. Therefore, it is essential to wait for the cooling time to pass before taking the parts out, otherwise there is a **risk of burns**.

Inadequate or ineffective disinfection encourages the growth of bacteria and thus increases the **risk of infection**. Disinfection is not completed properly unless both the specified minimum disinfection time and the required power output of the microwave are observed. The disinfector must also be perfectly clean. Check regularly that the microwave and the disinfector are functioning perfectly.

Use only water for disinfection, never use chemicals.

Chemical cleaning with disinfection

Cleaning and disinfection can be carried out in a single cycle using a chemical preparation process. For this, use the cleaning disinfecting agent Bomix® plus.

To ensure safety when handling chemicals, follow the instructions for use of the disinfecting agent, particularly the accompanying safety instructions.

Inadequate disinfection encourages the growth of bacteria and thus increases the **risk of infection**. Adequate cleaning with disinfection can only be assured if the specified mixing ratio and application time are adhered to, and if all individual parts are completely immersed in the solution for the entire application time. There must not be any air pockets or bubbles.

- Prepare a 2% solution of Bomix[®] plus by mixing 10 ml of the concentrate with 500 ml tap water.
- Place all the individual parts in the prepared solution and leave them to soak for 5 min.

- If the application period is exceeded significantly, the plastic parts may take on the smell of the disinfectant.
- Rinse off all parts thoroughly in running water (residues of the disinfectant can cause allergic reactions or irritations of the mucous membrane).
- You can remove excess water more quickly by shaking all parts.
- Dispose of the used solution (the diluted solution can be got rid of down the drain).

4.5 Visual inspection

Inspect all product components after each cleaning and disinfection. Replace any broken, misshapen or seriously discoloured parts.

4.6 Drying and storage

- Place all parts on a dry, clean and absorbent surface and let them dry completely.
- Wrap the nebuliser in a clean, lint-free cloth (e.g., a teatowel) and keep it in a dry, dust-free environment.

5 HYGIENIC RE-USE IN THE PROFESSIONAL ENVIRONMENT

5.1 General

NOTICE

Ensure that the system is processed for re-use only with the application of adequately validated device- and product-specific procedures, and that the validated parameters are maintained during each cycle. The equipment used must be maintained regularly.

Follow the schedule defined below for hygienic re-use:

- Single patient use:
 - Clean the nebuliser and the accessories used after every inhalation session.
 - Disinfect the nebuliser and accessories used at least once a day.
- Multiple patient use:
 Carry out the following steps before every change in patients:

- Clean, disinfect and sterilise the nebuliser and the accessories used.
- Clean and disinfect the connection tubing mechanically, or replace it. For this procedure, follow the instructions in the leaflet "Mechanical cleaning with disinfection of PARI connection tubes" (Doc. no. 041D0503), available on the internet at www.pari.de.

PARI nebulisers must be processed for re-use with the application of suitable, validated methods in such manner that the success of these methods can be monitored reliably and the safety and health of patients is not threatened.

⚠ CAUTION

Ensure that the components are allowed to dry sufficiently after each processing step. Condensation or residual wetness can present an increased **risk of bacterial growth**.

The following text describes the methods validated by PARI, with which plastic PARI parts can be processed effectively.

When processing for re-use, follow the instructions in the leaflet "Hygienic re-use in the professional environment" (Doc. no. 041D0498), available on the internet at www.pari.de.

5.2 Processing limits

The nebuliser can be sterilised up to 300 times and used for up to one year.

Do not expose the nebuliser or accessories used to a temperature above 137 °C.

The connection tubing can be processed for re-use about 50 times

5.3 Preparation

- · Detach the tubing from the nebuliser.
- Make sure that all medication residues are removed from the nebuliser.
- · Dismantle the nebuliser into its individual parts.
- Carefully pull the blue expiratory valve out of the slot in the mouthpiece. The valve must still be attached to the mouthpiece.

Rinse all components of the nebuliser under running tap water. This must be done especially thoroughly if cleaning and disinfection are not carried out straight away.

5.4 Cleaning and disinfection

Generally, a mechanical process (instrument dishwasher) should be used for cleaning and disinfection.

NOTICE

When selecting the cleaning agent and/or disinfectant, ensure that it is compatible with the materials used, **otherwise** damage to the product components cannot be ruled out.

Information when using masks with rubber bands:
The rubber band can only be disinfected chemically, not thermally as it might be damaged by the high temperatures used.

Mechanical cleaning with disinfection

FQUIPMENT

- Cleaning and disinfection device in conformance with DIN EN ISO 15883
- Cleaning agent

PROCEDURE

- Arrange the components so that they can be cleaned and disinfected most effectively.
- · Select a suitable cleaning programme and agent.

VALIDATION

The method has been validated in Europe using:

- G7836 CD cleaner-disinfector produced by Miele (programme Vario TD²)
- Standard baskets produced by Miele for nebulisers and accessories (without the connection tubing) or special baskets from Miele for nebulisers, accessories and the connection tubing
- Neodisher Mediclean forte 0.5% (alkaline cleaning agent) and Neodisher Z (citric acid for neutralisation) from Dr. Weigert

DRYING

Even if the cleaner-disinfector is equipped with a drying function, make sure that no residual moisture is left in the components

- If necessary, remove any residual moisture by shaking the parts and allow them to dry completely.
- Allow the compressor or the PARI CENTRAL to continue blowing air through the tubing until all of the condensate has been blown out.

²⁾ Cleaning: 50 minutes at max. 65 °C; disinfection: 5 minutes at max. 93 °C

Manual cleaning

As an alternative to mechanical cleaning and disinfection, the following manual method can be used:

EQUIPMENT

- Cleaning agent
- Brush
- Warm water

PROCEDURE

- · Place all the parts in the cleaning solution.
- · Clean the parts thoroughly with a brush.
- Thoroughly rinse the parts under hygienically tested water and leave them to dry completely on a clean, dry and absorbent surface.

VALIDATION

The method has been validated using:

- A pH-neutral enzymatic cleaner, Korsolex[®] Endo Cleaner 0.5% (Bode)
- Approx. 50 °C warm water
- Duration of cleaning procedure: approx. 5 minutes

Manual disinfection

For manual disinfection, select a suitable instrument disinfectant that has been proven to be effective and compatible with the cleaning agent used.

EQUIPMENT

- Disinfectant
- Water

PROCEDURE

- · Place all the parts in the disinfection solution.
- Allow the solution to take effect as indicated in the instructions for use of the disinfectant.
- Thoroughly rinse the parts under hygienically tested water and leave them to dry completely on a clean, dry and absorbent surface.

VALIDATION

The method has been validated using an aldehyde-containing disinfectant:

- Korsolex® Basic 4% (Bode)
- Treatment time: 15 minutes

Manual cleaning with disinfection


The following chemical process may also be used as an alternative to the mechanical process:

FQUIPMENT

- Cleaning disinfectant (a combination product for cleaning and disinfection in a single work step)
- Warm water
- Basin

PROCEDURE

Prepare the cleaning and disinfectant solution in the basin in accordance with the manufacturer's instructions. The quantity of solution required depends on the number of nebulisers to be cleaned.

- Place all the individual parts in the prepared solution and leave them to soak.
- · Clean the single parts thoroughly with a brush if necessary.
- Thoroughly rinse all parts under hygienically tested water and leave them to dry completely on a clean, dry and absorbent surface.
- · Dispose of the used solution.

VALIDATION

The method has been validated using:

- cleaning instrument disinfectant Bomix[®] plus 2% (Bode)
- Treatment time: 5 minutes
 Info: If the application period is exceeded significantly, the plastic parts may take on the smell of the disinfectant.

Inadequate disinfection encourages the growth of bacteria and thus increases the **risk of infection**. Adequate cleaning with disinfection can only be assured if the specified mixing ratio and application time are adhered to, and if all individual parts are completely immersed in the solution for the entire application time. There must not be any air pockets or bubbles.

A CAUTION

Disinfectant residues on the product can cause allergic reactions or irritation of the mucous membrane!

The diluted disinfectant solution can be disposed of down the drain.


5.5 Sterilisation

NOTICE

A validated sterilisation process can only be carried out on products that have been cleaned and disinfected.

Information on sterilising the child/adult mask soft:

Always use the mask stabiliser when sterilising these masks, because otherwise they may lose their shape under the effects of high temperatures. For this purpose, follow the instructions for use of the mask stabiliser as well.

FQUIPMENT

Steam steriliser – preferably with fractionated pre-vacuum – in accordance with DIN EN 285 or DIN EN 13060 (type B)

PROCEDURE

Pack all of the disassembled parts in a sterile barrier system in accordance with DIN EN 11607 (e.g., foil-paper packaging).

Maximum sterilisation temperature and holding time:

121 °C, at least 20 min. or

132 °C / 134 °C, at least 3 min.

VALIDATION

The method has been validated according to DIN EN ISO 17665-1 using:

- Euro-Selectomat 666 from MMM (fractionated pre-vacuum method)
- Varioklav 400 E from H+P Labortechnik (gravitation method)

5.6 Visual inspection

Inspect all product components after each cleaning, disinfection and sterilisation. Replace any broken, misshapen or seriously discoloured parts.

5.7 Storage

Store the nebuliser and the accessories used in a dry, dustfree location equipped with protection against contamination. It is recommended to use sterile packaging.

6 MISCELLANEOUS

6.1 Disposal

All product components can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.

6.2 Technical data


Operating gases	Air, oxygen
Minimum flow	3.0 l/min.
Corresponding minimum operating pressure	0.5 bar / 50 kPa
Maximum flow	6.0 l/min.
Corresponding maximum operating pressure	2.0 bar / 200 kPa
Minimum fill volume	3 ml
Maximum fill volume	8 ml

6.3 Aerosol characteristics

In accordance with DIN EN 13544-1 Appendix CC.

Determination of particle size

Parameter	
Measurement time	3 min
Fill volume	3 ml
Test solution	Sodium fluoride (2.5%)
Test instrument	Marple cascade impactor
minimum flow	3 l/min.
maximum flow	6 l/min.

Averages from each of two individual measurement series with three nebuliser-compressed air combinations

Determination of aerosol output

The aerosol output is the quantity of aerosol that the nebuliser system produces with a given fill volume.

Parameter	
Fill volume	2 ml
Test solution	Sodium fluoride (1%)
Test instrument	PARI COMPASS breathing simulator
Measurement time	Until the exiting aerosol stream becomes irregular, plus 1 minute

Aerosol output rate ³	
at minimum flow (3 l/min)	0.14 ml
at maximum flow (6 l/min)	0.15 ml

Determination of the aerosol output rate


The aerosol output rate is the quantity of aerosol that the nebuliser system produces in a given time.

Parameter	
Fill volume	2 ml
Test solution	Sodium fluoride (1%)
Test instrument	PARI COMPASS breathing simulator
Measurement time	1 minute

Aerosol output rate ³	
at minimum flow (3 l/min)	0.05 ml/min
at maximum flow (6 l/min)	0.11 ml/min

6.4 Explanation of symbols

The following symbols can be found on the nebuliser and/or the packaging:

6.5 Contact

For all product information and in the event of defects or questions about usage, please contact our Service Center:

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