
eBase[®] Controller

For Use with the eFlow[®] Technology
nebuliser handsets Tolero[®] for Vantobra[®]
and eFlow[®] rapid for Colfinair[®]

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
For full details, please read and take note of the instructions for use.

1 ASSEMBLE

* TouchSpray made under licence from The Technology Partnership plc.


Short Guide

For full details, please read and take note of the instructions for use.

2 INHALE


Sit in an upright position.
Hold nebuliser handset horizontally.


For full details, please read and take note of the instructions for use.

3 NEBULISER HANDSET-HYGIENE AT HOME


Clean immediately after inhalation

Warm water, clear liquid dish soap

Disinfect

Disinfect in the steam disinfectant for 6 minutes, or boil thoroughly in distilled water for 5 minutes.

Dry

Allow to air-dry completely. Store in a clean, dry place.

en	Instructions for use	3
el	Οδηγίες χρήσεως	31
es	Instrucciones de uso	61
it	Istruzioni per l'uso	89

For eBase® Controller Type 678G1020

Table of contents

1	IMPORTANT INFORMATION	5
	Organisation of warnings	5
	Handling the nebuliser system	5
	Treatment of babies, children and individuals requiring assistance	6
2	PRODUCT DESCRIPTION	7
	Intended purpose	8
	Indication	8
	Intended use	8
	Contraindications	8
	Intended patient groups	8
	Lifetime of the medical device	8
3	PREPARING FOR INHALATION	9
	Information for the first application	9
	Preparing the controller	9
	Mains power	9
	Preparing the nebuliser handset	10
4	PERFORMING INHALATION	13
	Pausing during inhalation	13
	Screens during inhalation	14
	End of the inhalation session	14
5	CLEANING AND DISINFECTION	16
	Nebuliser handset	16
	Preparation	16
	Hygienic reprocessing at home	17
	Cleaning	17
	Disinfection	18
	Hygienic reprocessing in the hospital and doctor's practice	19
	Cleaning with thermal disinfection	19
	Controller, connection cord and power adapter	20
	Cleaning	20
	Disinfection	20
	Storage	20

6 TROUBLESHOOTING 21

7 DISPOSAL 24

8 TECHNICAL DATA 24

 General 24

 Classification according to EN 60601-1 25

 Materials used 26

 Operating conditions 26

 Storage and transportation conditions 26

9 SYMBOLS 27

10 TERMS AND CONDITIONS OF WARRANTY 29

 Service and manufacturer 29

1 IMPORTANT INFORMATION

Read these instructions for use and the instructions for use of the accompanying accessories completely before using for the first time. Keep them safely so that you can refer to them again later.

CAUTION:

If these instructions for use are not followed, injuries or damage to the product cannot be ruled out.

Report serious incidents to the manufacturer and the competent authority.

You should have a medical examination before treating illnesses.

Organisation of warnings

In these instructions for use, safety-related warnings are organised in hazard levels:

- The signal word **WARNING** is used to indicate dangers that may result in serious injuries or even death if precautionary measures are not taken.
- The signal word **CAUTION** is used to indicate dangers which may result in slight to moderate injuries or impair treatment if precautionary measures are not taken.
- The signal word **NOTE** is used to indicate general precautionary measures for working with the product which should be followed to avoid damaging the product.

Handling the nebuliser system

Check the nebuliser system before every application.

CAUTION:

Replace all broken, deformed and heavily discoloured parts. Damaged parts can impair the function of the nebuliser system and thus also affect treatment.

WARNING:

Do not operate the nebuliser system if the power adapter is visibly damaged, since otherwise there is a danger from contact with live parts (e.g., electric shock).

Follow the instructions below in order to operate the nebuliser system safely:

- Always remove the power adapter from the socket to ensure that the power is completely cut off.

WARNING:

Never remove the power adapter from the socket with wet hands. There may be a danger of electric shock.

- Keep the nebuliser system away from hot surfaces (e.g., hob).
- Do not allow pets (e.g., rodents) near the cables.

WARNING:

Do not use the nebuliser system in areas where there is a risk of explosion or in the presence of gases that promote combustion (e.g., oxygen, nitrous oxide, flammable anaesthetics).

- Do not carry out an inhalation session while you are driving a motor vehicle (risk of accident).

Treatment of babies, children and individuals requiring assistance


WARNING:

Children and anyone who requires assistance must be supervised constantly by an adult during inhalation therapy. This is the only way to ensure safe and effective treatment. These individuals often underestimate the hazards associated with the equipment (e.g. risk of strangulation with the cable), and this can result in risk of injury.

WARNING:

The product contains small parts. Small parts can block the airways and pose a choking hazard. Therefore, be sure to keep them out of the reach of children at all times.

2 PRODUCT DESCRIPTION


- 1 eBase® Controller
- 1a ON/OFF button
- 2 Connection cord (connection between controller and nebuliser handset)
- 3 Power adapter
- 3a Country-specific adapter (AU, EU, UK, US)
- Carrying case
- Batteries (optional)

Check that all components are contained in your package. If they are not, contact your service partner.

Intended purpose

The eBase[®] Controller is used to operate eFlow[®] Technology nebuliser handsets, e.g., Tolero[®] for Vantobra[®] and eFlow[®] rapid for Colfinair[®]. The nebuliser handset is not included in this kit, and is made available for use with a specially approved drug and/or via a specific distribution channel.

Indication

The eBase[®] Controller has no indication.

Intended use

The eBase[®] Controller must only be operated in accordance with its intended purpose.

WARNING:

For reasons of hygiene and to avoid cross-infections, the nebuliser system (the eBase[®] Controller in conjunction with a nebuliser handset) must only be used by one patient.

WARNING:

Do not use any other drug than the one approved for the nebuliser handset.

CAUTION:

Never attempt to inhale essential oils, as this can cause irritation of the airways (intense coughing).

Only original PARI accessories must be used to operate the system.

CAUTION:

Using accessories and spare parts other than those supplied by the manufacturer may cause increased electromagnetic emissions or reduced electromagnetic immunity, leading to a malfunction.

PARI Pharma GmbH is not responsible for damage if the system is used improperly or contrary to the instructions.

Contraindications

The eBase[®] Controller has no contraindication.

Intended patient groups

The nebuliser system guarantees highly effective, fast inhalation treatment for patients aged from infant to adult. The exact age at which treatment can start depends on the drug and may vary from one drug to another.

The nebuliser system is only suitable for use by patients who can breathe unassisted and are conscious.

i *Info:*
Besides the information in these instructions for use, please also take note of the package inserts for the nebuliser handset and the medication.

Lifetime of the medical device

- The expected life time of the eBase[®] Controller is 3 years.
- Nebuliser handset: Follow the instructions for use of the nebuliser handset you are using.

3 PREPARING FOR INHALATION

Information for the first application

- Wash your hands thoroughly every time before using the nebuliser system.
- To avoid irritation of the airways due to inhaling aerosol when it is too cold and impairment of the aerosol properties, do not use medications straight from the refrigerator.

WARNING:


Clean the nebuliser handset and aerosol head before the first application (see "CLEANING AND DISINFECTION", page 16).

Preparing the controller


The controller can be operated with batteries or with rechargeable batteries.

i *Info:*
Rechargeable batteries cannot be recharged in the controller!

- Open the battery compartment by lifting the clip.

- Insert the batteries as indicated by the polarity symbols.

- Close the battery compartment.

Mains power


Check that the power adapter is intact before every application.

WARNING:


Do not use a power adapter that is damaged or faulty.

The eBase[®] Controller was developed for mobile use. Therefore, insert batteries before use and leave them in the controller during mains operation. This guarantees smooth functionality. This does not deplete the energy reserves of the batteries.

- Replace the country-specific adapter (3a) for the power adapter as required. To unlock, press the clip and push the country-specific adapter upwards. When plugging in, make sure that the adapter engages in the power adapter housing.

- Connect the power adapter to the power cord socket on the controller.

When the adapter is plugged into the socket, the controller is ready for use.

Preparing the nebuliser handset



CAUTION:

Check the parts of your nebuliser handset before every application. Replace all broken, deformed and heavily discoloured parts. Please also follow the instructions for assembly given below. Damaged parts and incorrectly assembled nebuliser handsets can impair the function of the nebuliser handset and thus also affect the treatment.


- Insert the aerosol head (4) into the openings on the medication reservoir (5) (opened aerosol chamber (5a)).

CAUTION:


Do not touch the membrane (the curved metal surface in the middle of the aerosol head), it can be damaged very easily.

- Insert the inspiratory valve (6). When doing so, check that the two valve flaps are seated correctly. They must lie perfectly flush and face towards the aerosol chamber following insertion.


- Close the aerosol chamber.


NOTE:

If the snap closure does not close, ensure that the inspiratory valve is seated correctly. It must lie flat against the aerosol chamber.


- Attach the mouthpiece (7) to the aerosol chamber.

- Connect the connection cord (2) to the nebuliser handset.


- Connect the connection cord (8) to the controller.


- Add the quantity of medication prescribed by your doctor to the medication reservoir.

WARNING:

Ensure that the medication is filled no higher than the upper scale graduation mark (6 ml). Do not overfill under any circumstances!


- Check that all parts are connected tightly and that the medication reservoir is sealed.

i *Info:*
If you intend to inhale several medications one after the other, rinse all nebuliser handset parts thoroughly with warm tap water before refilling.

- Close the medication reservoir.

4 PERFORMING INHALATION

- Sit in an upright position and relax. This makes inhaling easier and improves deposition of the medication in the airways.
- Hold the nebuliser handset in your hand.


CAUTION:

Hold the nebuliser handset horizontally throughout the inhalation session. If held in a tilted position, an incorrect dosage may occur.

- Hold the mouthpiece between your teeth and enclose with your lips. The lips should not touch the blue expiratory valve.

- Press the ON/OFF button on the controller to start the inhalation session.

A green LED next to the ON/OFF button lights up and an audible signal (1 beep) is emitted to indicate readiness for use.

- Inhale and exhale slowly and deeply through the mouthpiece. The mouthpiece should also remain in the mouth while exhaling. Do not breathe through your nose. You may use a nose clip if your doctor agrees.

i Info:
It is normal for aerosol to escape through the expiratory valve in the mouthpiece during exhalation; this is not a malfunction.

Pausing during inhalation

If it is necessary to interrupt the inhalation session briefly, use the pause mode:

- To activate the pause mode, press the ON/OFF button for about 1 second. Pause mode is activated if the LED flashes green.
- To continue inhalation, press the ON/OFF button for about 1 second once again.

i Info:
Pause mode can only be activated after an operating time of a few seconds.

If pausing for a longer time, switch the controller off by pressing the ON/OFF button for approx. 2 seconds.

As long as there is enough liquid in the medication reservoir, the nebuliser handset can be restarted by pressing the ON/OFF button again.


Screens during inhalation

The following screens appear during the inhalation session:


Start

During inhalation


Nebulising with battery power


Nebulising with mains power

Activated pause mode


Pause mode during battery operation

Pause mode during mains operation


End of the inhalation session

The inhalation session is finished as soon as no more aerosol is generated, and the following screen appears:

The controller is switched off automatically. If the controller switches itself off considerably sooner or later than usual, follow the instructions in Chapter "TROUBLESHOOTING", page 21.

i Info:
Note that with the eFlow[®]rapid nebuliser handset about 1 ml of medication cannot be nebulised and remains unused in the medication reservoir; this residue must be discarded. This is necessary to ensure a safe quantity of medication is delivered and does not constitute a malfunction.

For safe disposal of the leftover medication, follow the information in the medication package leaflet. The nebulisation time depends on the quantity of medication added, and varies according to the medication used.

 **Info:**

I *If the nebuliser system is used regularly, the nebulisation time may slowly increase. This is not a malfunction.*

- When operating with the power adapter: When the inhalation session has ended, disconnect the power adapter plug from the socket.
- Carry out the hygiene procedures (see "CLEANING AND DISINFECTION" on page 16) immediately after the inhalation session has ended.

5 CLEANING AND DISINFECTION

Nebuliser handset

WARNING:

The nebuliser handset (and the aerosol head) must be cleaned immediately after each use and disinfected once a day.

Preparation

Dismantle the nebuliser handset into its individual parts:


- Detach the nebuliser handset from the connection cord adapter.
- Remove the mouthpiece from the nebuliser handset.
- Open the medication reservoir and pour out the medication residue.

i Info:

The medication reservoir and aerosol chamber are already preassembled and must not be detached from one another.


- Open the aerosol chamber (5a).

- Remove the inspiratory valve from the aerosol chamber.
- Remove the aerosol head (4). Press lightly on the engaging hooks on the sides of the aerosol head.

CAUTION:

Do not touch the membrane (the curved metal surface in the middle of the aerosol head) to avoid damage.

Hygienic reprocessing at home

Cleaning

WARNING:

Clean the nebuliser handset (incl. aerosol head) immediately after each application.


NOTE:

Rinsing all the parts of the nebuliser handset and the aerosol head in warm tap water straight after use prevents medication residue from accumulating, for example. This helps you to prevent potential extension of the nebulisation time due to partial closure of the pores in the membrane.


- Rinse each side of the aerosol head for 5 seconds each in warm tap water (water of drinking water quality, temperature up to 37 °C).
- Rinse all parts of the nebuliser handset for approx. 5 seconds.
- Place all nebuliser handset parts as well as the aerosol head in warm tap water with a few drops of clear liquid dish soap for 5 minutes. Use a clean bowl for this.

CAUTION:

Do not use a brushes or any other tools to clean the aerosol head as they may cause damage to the product components.

- Rinse all nebuliser handset parts and the front and back of the aerosol head thoroughly under running tap water.
- Shake all parts to remove excess water more quickly.
- Place all nebuliser handset parts on a dry, clean surface and allow to dry completely.

CAUTION:

Moisture promotes the growth of bacteria. Please ensure adequate drying after each cleaning.

Disinfection

WARNING:

Disinfect the nebuliser handset parts and the aerosol head at least once a day, e.g., immediately after cleaning at the end of the day (a nebuliser handset cannot be disinfected effectively unless it has been cleaned first).

NOTE:


Do not use a microwave to disinfect the aerosol head. The possibility of product components being damaged cannot be ruled out.

Using a standard commercial disinfecter for baby bottles (so-called steam sterilizer)

- To ensure effective disinfection, use a thermal disinfecter with an operating time of at least 6 minutes. With regard to performing the disinfection, the duration of the disinfection process and the quantity of water required, follow the instructions for use of the disinfecter you are using.

CAUTION:

Inadequate disinfection promotes the growth of bacteria and thus increases the risk of infection. If disinfection is not carried out for the time specified for the device concerned, it will not be effective. Therefore, do not switch the device off before the disinfection cycle is complete. Also, make sure that the device is clean and check it regularly for proper operation.

- Take the parts out of the disinfecter immediately after disinfection is completed. Place them on a dry, clean surface and allow them to dry completely.

CAUTION:

Check the nebuliser handset parts after every disinfection. Replace all broken, deformed and heavily discoloured parts.


Possible option: Boiling in distilled water

- Place the individual parts of your nebuliser handset (incl. aerosol head) in boiling water for at least 5 minutes. Use a clean pot and distilled water for this.

CAUTION:

Plastic will melt if it comes in contact with the hot bottom of the pot. Therefore, make sure that there is sufficient water in the pot. This will prevent the parts from being damaged.

- Take the parts out of the pot.
- Shake all parts to remove excess water more quickly.
- Place all nebuliser handset parts on a dry, clean and surface and allow to dry completely.

CAUTION:

Moisture promotes the growth of bacteria. Please ensure adequate drying after every disinfection.

CAUTION:

Check the nebuliser handset parts after every disinfection. Replace all broken, deformed and heavily discoloured parts.

Hygienic reprocessing in the hospital and doctor's practice

Carry out cleaning with disinfection immediately after use. Ideally, a mechanical procedure (washer-disinfector) should be used.

Prepare for cleaning with disinfection (see "CLEANING AND DISINFECTION" on page 16).

Cleaning with thermal disinfection

- Place the nebuliser handset components in the washer-disinfector.
- Select the 93 °C programme (treatment time: 10 minutes).

The effectiveness of this procedure has been validated using a G7736 disinfector made by Miele and using neodisher MA (Chemische Fabrik Dr. Weigert, Hamburg) as a cleaning agent in conjunction with deionised water as a neutraliser.

Even when the washer-disinfector is equipped with a drying function, ensure that no residual moisture remains in the nebuliser handset. Remove residual moisture by shaking the parts and allow them to air-dry completely.

Controller, connection cord and power adapter

Cleaning

Clean the controller, the connection cord and the power adapter as necessary with a clean, lint-free, damp cloth.

Allow the controller, the connection cord and the power adapter to dry completely.

NOTE:

Do not hold the controller, the connection cord or the power adapter under running water and do not use liquid cleaning agents.

If liquid gets into these components, the electronics may be damaged and the device may no longer function properly.

If liquid does get into them, contact the local service partner immediately.

Disinfection

Disinfect the controller, the connection cord and the power adapter after cleaning. Use an alcohol-based, standard commercial disinfectant wipe for disinfecting. To ensure safe application, follow the instructions for use of the disinfectant wipe. Effectiveness has been tested using Bacillol Tissues and Clinell Wipes.

Allow the controller, the connection cord and the power adapter to dry completely.

Storage

Between uses, and in particular during long breaks in treatment, keep the nebuliser handset in a dry, dust-free environment (e.g., nebuliser handset bag).





Pack the controller, the power adapter and the connection cord in the carrying case provided for this purpose.

NOTE:

Leaking batteries may cause damage to the controller. Therefore, remove the batteries or rechargeable batteries if the nebuliser system is not expected to be used for a relatively long period of time.

6 TROUBLESHOOTING

You can correct malfunctions that occur during daily operation yourself. The following tables will help you to locate and correct the causes of malfunction.

Screen (flashes alternately)	Possible cause	Solution
	Connection between the controller and the nebuliser handset/aerosol head interrupted.	<ul style="list-style-type: none"> - Verify that all parts are securely connected. - Replace the aerosol head and/or the connection cord if damaged. If the fault persists, contact your service partner.
	No medication added.	Add medication.
	Inhalation interrupted.	Resume inhalation by restarting the device (press ON/OFF button).
	The maximum operating time of 20 minutes for an inhalation session has been exceeded.	Resume inhalation by restarting the device (press ON/OFF button).
	Battery charge level 50%.	Have new batteries on hand.
	Battery flat.	Insert new batteries or operate with the power adapter connected.

Fault/Condition	Possible cause/Solution
<p>The controller cannot be activated (no beep, no orange or green light signal).</p>	<ul style="list-style-type: none"> - To switch on, press the ON/OFF button for approx. 2 seconds. - Is power cord plug plugged into mains and the device (for mains operation)? - Are batteries or rechargeable batteries inserted correctly? - Check the charge level of the batteries.
<p>After switching the controller on, no aerosol is generated, and/or the controller switches off again after a few seconds.</p>	<ul style="list-style-type: none"> - Has the medication reservoir been filled with medication? - Check the connection between the controller and the nebuliser handset.
<p>The controller switches off, even though there is still medication in the medication reservoir.</p> <p>i <i>Info:</i> <i>Note for the eFlow[®]rapid nebuliser handset that about 1 ml of medication cannot be nebulised and remains unused in the medication reservoir; this residue must be discarded. This is intentional and does not represent a malfunction. If the Tolero[®] nebuliser handset is used no significant residual volume remains in the medication reservoir.</i></p>	<p>If there is still more than 1 ml in the medication reservoir of the eFlow[®]rapid nebuliser handset:</p> <ul style="list-style-type: none"> - Check the charge level of the batteries or rechargeable batteries! - The nebuliser handset was not being held horizontally. Hold the nebuliser handset horizontally. - The controller switches off automatically after 20 minutes. <p>Carry out inhalation to the end by pressing the ON/OFF button again, so that you can inhale the prescribed dose of medication.</p> <p>If the maximum operating time for each inhalation session is exceeded repeatedly, the aerosol head must be replaced.</p>
<p>The controller does not switch off automatically even when no more aerosol is being produced and there is no more than about 1 ml medication left in the medication reservoir if a eFlow[®]rapid nebuliser handset is used.</p>	<p>The treatment session has been completed successfully and the controller can be switched off by pressing the ON/OFF button.</p>

Fault/Condition	Possible cause/Solution
The nebulisation time has become noticeably longer.	<p>If the inhalation time with the same fill volume and the same substance becomes considerably longer.</p> <ul style="list-style-type: none"> - You are inhaling a different medication. - The fill volume has been increased. - The medication cap was not closed tightly. - The aerosol head was not cleaned immediately after the last application (the fine pores of the membrane may be partially blocked), or it is damaged. - You are using the aerosol head for longer than its recommended lifetime.
Aerosol escapes continuously in large volumes from the opening slots in the medication reservoir.	<ul style="list-style-type: none"> • Check whether the nebuliser handset has been assembled correctly. When doing so, check that the two valve flaps of the inspiratory valve are seated correctly. They must lie flat and face towards the aerosol chamber.
The medication cap cannot be removed from the medication reservoir after the inhalation treatment.	<p>The negative pressure in the medication reservoir is too great.</p> <ul style="list-style-type: none"> • Open the aerosol chamber and gently pull the aerosol head so that air is able to flow in. <p>The medication cap can now be removed easily.</p>

If the malfunction persists, contact your service partner.

7 DISPOSAL

Live parts (such as the controller, the cable, the power adapter and the batteries) and the aerosol head must not be

disposed of with normal domestic waste. The device components must be disposed of in accordance with local regulations.

8 TECHNICAL DATA

General

Electrical connection	PARI Power power adapter REF 078B7114 Input: 100 to 240 V~, 50-60 Hz Output: 12 V \equiv
Batteries	4 x 1.5 V Mignon (AA LR6/Alkaline)
Rechargeable batteries	4 x 1.2 V Mignon (AA rechargeable, at least 2100 mAh)
Weight of controller (incl. batteries)	approx. 300 g
Controller housing dimensions	H 4.0 cm, Ø 11.6 cm

Technical data on electromagnetic compatibility are available in table form upon request from PARI Pharma GmbH, or on the Internet at <https://www.pari.com/fileadmin/Electromagnetic-compatibility-4.pdf>

Classification according to EN 60601-1

Type of electric shock protection (power adapter)	Protection class II
Degree of protection from electric shock of applied part (nebuliser handset)	Type BF
Degree of protection against ingress of water as per IEC 60529 (IP rating)	IP 21
Degree of protection when used in the presence of flammable mixtures of anaesthetics with air or with oxygen or nitrous oxide	No protection
Operating mode	Continuous operation
Environment	<ul style="list-style-type: none"> - Professional healthcare facility environments - Domestic healthcare settings

Materials used

This product consists mainly of the following materials:

- Controller:
thermoplastic elastomers, acetonitrile butadiene styrene, polypropylene
- Nebuliser handset (incl. aerosol head):
polypropylene, thermoplastic elastomers, polyoxymethylene, silicone rubber, stainless steel
- Mouthpiece:
polypropylene, silicone rubber

The nebuliser handset contains no parts made from natural rubber (latex).

Operating conditions

- Ambient temperature:
+5 °C to +40 °C
- Relative humidity of ambient air:
15% to 93% (non-condensing)
- Air pressure:
700 hPa to 1060 hPa

Users can operate the device both at home and when out and about. Its use on trains and aeroplanes is restricted to passenger areas only. If using the device in a car, it must be operated with battery power.

Within professional healthcare facility environment, the device must only be used in in-patient wards and intensive care units. The device must never be used in areas exposed to elevated magnetic or electrical radiation (e.g. in the vicinity of an MRI scanner).

CAUTION:

When using portable wireless communication devices, you must ensure that a distance of at least 30 cm is maintained between all parts of the nebuliser system and these devices (including their accessories, e.g. antenna cables or external antennas). Otherwise, the performance of your nebuliser system may deteriorate.

CAUTION:

The controller must not be placed directly beside or on top of other devices for operation. If operation beside or on top of other devices is necessary, the controller must be monitored constantly to ensure proper operation.

Storage and transportation conditions

- Temperature: -25 °C to +70 °C
- Relative humidity of ambient air:
0% to 93% (non-condensing)
- Air pressure: 500 hPa to 1060 hPa

NOTE:

Moisture from condensation can impair the controller's functional capabilities. Therefore, you should avoid extreme temperature fluctuations. Wait until the controller has reached room temperature before beginning an inhalation session.

If the device is not expected to be used for a prolonged period, the batteries or rechargeable batteries must be removed.

Info:

Batteries vary considerably in terms of storage life and output depending on battery quality. The lifetime of about 90 minutes can only be achieved with the delivered batteries (4x).

9 SYMBOLS

	Medical Device
	Follow the instructions for use
	Single patient – multiple use
	Direct current
	Alternating current
	Lot number
	Item number
	Serial number
	Unique Device Identifier
	Protection class of the applied part: Type BF
	Protection class II equipment
IP21	The device is protected against dripping water
	Minimum and maximum ambient temperature
	Minimum and maximum humidity
	Minimum and maximum air pressure
	CE marking: This product complies with the European requirements for medical devices (EU 2017/745).
	The medical device was distributed commercially after 13 August 2005. This product must not be disposed of with normal domestic waste. The symbol of a crossed-through refuse bin indicates that separate collection is required.

Caution

Manufacturer

Date of manufacture

10 TERMS AND CONDITIONS OF WARRANTY

PARI guarantees to you that if used for its intended use your nebuliser system will be free from defects in material or workmanship caused by the manufacturing process for the warranty period of two years.

However this warranty does not cover wearing parts, that is to say device parts that are exposed to normal wear, such as the nebuliser handset and aerosol head.

The warranty shall no longer be valid if:

- the device has been operated or used in a manner inconsistent with the guidelines in the instructions for use,
- damage is evident that is attributable to external factors such as water, fire, lightning strike or similar,
- the damage has been caused by transporting the device incorrectly, or resulted from a fall,
- the device has been handled or serviced incorrectly,

- the serial number on the device has been changed, removed or rendered illegible,
- repairs, adaptations or modifications have been made to the device by persons not authorised by PARI.

If, exceptionally, a defect does occur, PARI will replace the device. The device may be replaced with an identically or at least comparably equipped model.

Replacement of the device shall not serve as grounds for new claims under the warranty. All exchanged used devices or parts will become the property of PARI. All other claims under the warranty are excluded to the extent permitted by law.

In the event of a claim under the warranty, please contact your local service partner. The original purchase receipt from the dealer must be presented as proof of warranty and ownership. This warranty applies solely for the initial purchaser of the device.

Service and manufacturer

Manufacturer	PARI Pharma GmbH Moosstr. 3 82319 Starnberg Germany
Contact for technical problems or questions regarding the device	Local service partner Contact: www.pari.com

The device must only be serviced by PARI Pharma GmbH or a service location expressly authorised by PARI Pharma GmbH.

We reserve the right to make technical changes.

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