


PEP® System I und II


CE 0123

de	Gebrauchsanweisung	1
en	Instructions for use	19
da	Brugsanvisning.....	37
el	Οδηγίες χρήσεως	55
es	Instrucciones de uso.....	73
fi	Käyttöohje	91
fr	Notice d'utilisation.....	109
it	Istruzioni per l'uso.....	127
nl	Gebruiksaanwijzing	145
no	Bruksanvisning.....	163
pt	Manual de Instruções	181
sv	Bruksanvisning.....	199
tr	Kullanım talimatları.....	217

A**B**

C

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Stand der Information: November 2015

en Instructions for use

Read these instructions carefully and completely before using the product for the first time. Keep these instructions in a safe place. Personal safety cannot be guaranteed if these instructions are not strictly followed.

Table of contents

1	Components	20
2	Intended use	21
3	Information	21
4	Preparation for treatment	24
4.1	Assembling the PARI PEP® System I	24
4.2	Assembling the PARI PEP® System II	25
4.3	Adjusting expiratory resistance	26
5	Application	27
5.1	Treatment without PARI nebuliser	27
5.2	Treatment in combination with PARI nebuliser	28
6	Hygienic re-use	29
6.1	General	29
6.2	Preparation	30
6.3	At home	31
6.4	Hospital and medical practice	32
7	Miscellaneous	35
7.1	Material resistance	35
7.2	Storage and transportation conditions	35
8	Explanation of symbols	36

1 Components

Figure A:

Please check the following table to ensure that all components listed are included in the product package. If a part is missing, contact the dealer from which you purchased the PARI PEP System immediately.

PARI PEP System I

(1) Mouthpiece without expiratory valve	(4) PEP II lower part ^{b)}
(2) PEP upper part	(5) Nose clip
(3) PEP I lower part ^{a)}	

PARI PEP System II

(1) Mouthpiece without expiratory valve	(6) Pressure gauge
(2) PEP upper part	(7) Adapter
(3) PEP I lower part ^{a)}	(8) PEP tubing adapter
(4) PEP II lower part ^{b)}	(9) Connection tube for pressure gauge
(5) Nose clip	

a) without pressure gauge attachment point

b) with pressure gauge attachment point

2 Intended use

The PARI PEP System provides an effective combination of PEP (Positive Expiratory Pressure) and inhalation treatment for secretion clearance in acute and chronic diseases of the lower respiratory tract. PEP treatment may be performed separately from or together with inhalation treatment.

The PARI PEP system can be used with **all** PARI nebulisers that are currently available commercially (with the exception of the nebulisers in the PARI LC SPRINT family, the PARI BABY and PARI TIA nebulisers). The manufacturer will be able to tell you whether the PARI PEP System can be used with PARI nebulisers that are introduced in the future.

Please read and follow the instructions for use of the corresponding PARI compressor and/or PARI nebuliser.

3 Information

The PARI PEP System is available in two versions, the PARI PEP System I without pressure gauge and the PARI PEP System II with pressure gauge. The gauge is used to adjust the expiratory resistance to your individual needs and to monitor it during treatment. Treatment with the PARI PEP System normally leads to an increase in lung volume and expiratory pressure, which in turn causes the interbronchial collateral channels to dilate. This promotes secretion clearance and, together with inhalation, ensures that more medication is delivered to the treatment site. The combination of PEP and inhalation treatment serves to clear secretions more quickly and efficiently. It also helps with exercising the expiratory and inspiratory muscles.

- You must always have a medical examination before treating illnesses. The PARI PEP System must not be used until you have been shown how to use it by your physician or physiotherapist.
- If you experience dizziness, shortness of breath or headaches during PEP therapy, stop the session immediately and notify your physician or physiotherapist.
- You must follow the general hygiene steps for using your PARI PEP System strictly, e.g. washing your hands and disinfection if required.

Note on the (optional) pressure gauge:

- The pressure gauge can be purchased as an optional accessory. It may be damaged if it is dropped or receives a similar impact. As a result, the expiratory resistance reading may be incorrect. If this occurs, it should be returned to PARI GmbH where it will be tested for a fee.
- The pressure gauge is used to test expiratory resistance, and as such has a measuring function. Accordingly, it is recommended that the device be calibrated once every 3 years. Please contact PARI GmbH or an authorised Service Centre for this.
- To prevent any risk to health from possible bacterial contamination of the pressure gauge connection tube (9), always use the PEP tubing adapter (8).

Material information

Polypropylene	Upper part, lower part, mouthpiece
Silicone	Upper part, tubing adapter
Thermoplastic elastomer	Nose clip
Polyacetal	Nose clip

Contraindications and side effects

The PARI PEP System must not be used by individuals who are suffering from untreated pneumothorax or actively coughing up blood.

The following side effects can occur:

dizziness, shortness of breath, headaches.

If you notice any of these side effects, discontinue the treatment immediately.

Contact your doctor or therapist to discuss how to proceed with treatment.

Inform your doctor or therapist if you notice any side effects that are not listed in these instructions for use.

4 Preparation for treatment

- Please follow the hygiene regulations before the device is used for the first time, before it is re-used, and after any long break in use.
- Ensure that the PARI PEP System was cleaned, disinfected and, if appropriate, sterilised after the last treatment (see Chapter 6 "Hygienic re-use").

4.1 Assembling the PARI PEP® System I

4.1.1 Treatment without PARI nebuliser

Figure A:

Follow these steps to assemble the PARI PEP System:

- Place the PEP upper part (2) on the PEP I lower part (3).
- Place the mouthpiece without expiratory valve (1) over the other end of the PEP I lower part.

4.1.2 Treatment in combination with PARI nebuliser

Figure C:

Follow these steps to assemble the PARI PEP System:

- Remove the inspiratory valve (10a) and the mouthpiece with the expiratory valve (10b) from the PARI nebuliser (10).
- Place the PEP upper part (2) on the PEP I lower part (3) and then place both together on the PARI nebuliser.
- Place the included mouthpiece without expiratory valve (1) over the PARI nebuliser.

4.2 Assembling the PARI PEP® System II

4.2.1 Treatment without PARI nebuliser

Figure A:

Follow these steps to assemble the PARI PEP System:

- Place the PEP upper part (2) on the PEP II lower part (4).
- Place the mouthpiece without expiratory valve (1) over the other end of the PEP II lower part.
- Place the PEP tubing adapter (8) on the pressure gauge attachment point of the PEP II lower part.
- Connect the connection tube of the pressure gauge (9) using the adapter (7) with the PEP tubing adapter.
- Insert the connection tube on the pressure gauge (6).

4.2.2 Treatment in combination with PARI nebuliser

Figure C:

Follow these steps to assemble the PARI PEP System:

- Remove the inspiratory valve (10a) and the mouthpiece with the expiratory valve (10b) from the PARI nebuliser (10).
- Place the PEP upper part (2) on the PEP II lower part (4) and then place both together on the PARI nebuliser.
- Place the included mouthpiece without expiratory valve (1) over the PARI nebuliser.
- Place the PEP tubing adapter (8) on the pressure gauge attachment point of the PEP II lower part.
- Connect the connection tube of the pressure gauge (9) using the adapter (7) with the PEP tubing adapter.
- Insert the connection tube on the pressure gauge (6).

4.3 Adjusting expiratory resistance

The pressure gauge (6) can be used to adjust expiratory resistance and to monitor the resistance constantly throughout PEP therapy. You must follow the instructions of your physician or physiotherapist for adjusting the expiratory resistance.

- Follow the steps described in Chapter 4.2 to assemble the PARI PEP System.

To obtain the best possible adjustment of expiratory resistance, you must not breathe through your nose. To prevent this, please use the nose clip (5).


Figure B:

Adjust expiratory resistance as follows, as instructed by your physician or physiotherapist:

- Remove the PEP upper part (2) until the engaging notch "x" releases catch "y" on the PEP II lower part (4).
The PEP upper part can now be easily twisted in both directions.
- Line up the appropriate opening in the PEP upper part opposite the opening in the PEP lower part.
- Push the PEP upper part onto the PEP lower part until the engaging notch "x" for the selected aperture is engaged with catch "y".
- Close your lips firmly around the mouthpiece (1) and breathe in. When you breathe out, you can read the expiratory resistance on the pressure gauge.

If the pressure displayed is not the pressure you want, choose a smaller or larger aperture as directed by your physician.

- Repeat this procedure until you achieve the desired expiratory resistance. Make a note of the corresponding aperture (see the markings on the PEP upper part 1.5; 2.0; 2.5; 3.0; 3.5; 4.0; 4.5 or 5.0).

5 Application


- Children and anyone who requires assistance must be supervised at all times by an adult during the PEP treatment session. Otherwise, safe and effective treatment cannot be guaranteed.
 - Cough to clear any secretion that is loosened during the treatment. Avoid coughing the secretions into the PARI PEP System.
- Caution:** Contains small parts which can be swallowed!

5.1 Treatment without PARI nebuliser


Refer to Chapter 4.1.1 or 4.2.1 to ensure that you have assembled the PARI PEP System correctly. Check that you have chosen the correct expiratory resistance setting (see Chapter 4.3). Use the nose clip to help you obtain the most efficient PEP treatment possible.

Perform the therapy as follows:

- Sit in an upright position and relax.
- Hold the mouthpiece (1) between your teeth, close your lips tightly around it and inhale.

- Exhale through the mouthpiece. The exhaled air should pass through the openings on the PEP upper part (2).

5.2 Treatment in combination with PARI nebuliser

Refer to Chapter 4.1.2 or 4.2.2 to ensure that you have assembled the PARI PEP System and attached the accompanying mouthpiece without expiratory valve correctly. Use the nose clip to help you obtain the most efficient PEP treatment possible. Please follow the instructions for use of your PARI nebuliser during an inhalation treatment session.

If there is still some medication left in the nebuliser at the end of the PEP treatment session, continue the inhalation treatment without the PARI PEP System. Note that it is normal for as much as 1 ml of medication to be left in the nebuliser.

6 Hygienic re-use

6.1 General

In order to prevent a risk to health (e.g. an infection occurring as a result of a contaminated PEP System) it is imperative that the following hygiene regulations be followed.

The PARI PEP System is intended for multiple use. Please be aware that different requirements regarding hygienic re-use apply to different fields of application:

- At home (use by a single patient):
The PARI PEP System, and where applicable the PEP tubing adapter, must be cleaned and disinfected after each use.
- Hospital/medical practice (possible use by more than one patient):
The PARI PEP System, and where applicable, the PEP tubing adapter, must be cleaned, disinfected and sterilised after each use.

The cleaning, disinfection, and sterilisation procedures described in Chapter 6.3 or 6.4 apply to all individual parts of the PARI PEP System except the pressure gauge, the pressure gauge connection tube and the adapter. These parts are to be wiped down with a damp cloth when necessary.

Additional requirements concerning the necessary hygienic preparations (hand care, handling of medications or inhalation solutions) in high-risk groups (e.g., patients with cystic fibrosis) can be obtained from the relevant support groups. Proof of suitability of the PARI PEP System for effective cleaning, disinfection and sterilisation has been obtained using an independent testing laboratory and using the recommended procedures stated in Chapter 6.4.1 "Cleaning/disinfection" (test documentation available for inspection on request). Using the specified alternatives remains the responsibility of the user.

Please ensure adequate drying after each cleaning, disinfection and/or sterilisation. Condensation or residual wetness can present an increased risk through the growth of bacteria.

Please check the parts of your PARI PEP System regularly and replace any defective (broken, misshapen, discoloured) parts.

Additionally, please read and follow the hygiene instructions in the instructions for use of the PARI compressor and the PARI nebuliser.


6.2 Preparation

Immediately after each treatment, the PARI PEP System must be cleaned of medication residues and contaminants. To do this, take the PARI PEP System completely apart (as described below).

Treatment without PARI nebuliser

Figure A:

- If necessary, disconnect the PEP tubing adapter (8) from the pressure gauge connection tube (9), and then disconnect it from the PEP II lower part (4).
- Pull the mouthpiece (1) and the PEP upper part (2) off the PEP lower part (3 or 4).
- Remove the valve shim (2a) from the PEP upper part.

Treatment in combination with PARI nebuliser

Figure C:

- If necessary, disconnect the PEP tubing adapter (8) from the pressure gauge connection tube (9), and then disconnect it from the PEP II lower part (4).
- Detach the PEP upper part (2) together with the PEP lower part (3 or 4) from the nebuliser.
- Detach the PEP lower part from the PEP upper part.
- Remove the valve shim (2a) from the PEP upper part (see figure above).

For hygienic re-use at home, please continue reading the following section.
For hygienic re-use in the hospital or medical practice, please read
Chapter 6.4 "Hospital and medical practice".

6.3 At home

The cleaning and disinfection procedure described in this chapter applies to all individual parts of the PARI PEP System with the exception of the pressure gauge, the pressure gauge connection tube and the adapter. These parts are to be wiped down with a damp cloth when necessary. Replace the tubing if it is contaminated.

6.3.1 Cleaning

- Take the PARI PEP System apart completely (see Chapter 6.2 "Preparation").

Recommended method: Use warm tap water

- Clean all parts thoroughly for 5 minutes with warm tap water (approx. 40 °C) and a little dishwashing liquid (in the strength specified in the manufacturer's recommendations).
- Then rinse all parts thoroughly under warm running water (approx. 40 °C without dishwashing liquid).
- Pour out the collected water.

You can remove excess water more quickly by shaking the parts.

6.3.2 Disinfection

Disinfect the dismantled PARI PEP System after cleaning.

In boiling water

- Place the individual parts in boiling water for at least 5 minutes. Use a clean pot and fresh drinking water.

If your water is very hard, limescale may be deposited on the parts as a milky coating. This can be prevented by using softened water.

Plastic will melt if it comes into contact with the hot base of the pot. Therefore, make sure there is plenty of water in the pot. This way, you will avoid **damaging the components**.

Using a standard thermal disinfector for baby bottles (not a microwave oven)

Disinfection must last at least 15 minutes. Regarding disinfection and the quantity of water required for this, it is essential to follow the instructions for use of the disinfector you are using. Always check carefully that it is clean and fully functional.

A damp environment may encourage the growth of bacteria. Therefore, remove the individual components from the pot or disinfecto as soon as disinfection has finished. Dry the parts as described in chapter 6.3.3. The **risk of infection** is reduced when the parts are dried completely.

6.3.3 Drying

- Place the components of the PARI PEP System on a dry, clean and absorbent surface and let them dry completely (at least 4 hours).

6.3.4 Storage

- Between use, and in particular during long breaks in treatment, wrap the PARI PEP System in a clean, lint-free cloth (e.g., a tea-towel) and keep it in a dry, dust-free environment.

Please continue reading from Chapter 7 "Miscellaneous".

6.4 Hospital and medical practice

The cleaning, disinfection, and sterilisation procedure described in this chapter applies to all individual parts of the PARI PEP System with the exception of the pressure gauge, the pressure gauge connection tube and the adapter. These parts are to be wiped down with a damp cloth when necessary. Replace the tubing if it is contaminated.

When the PARI PEP System is used for several patients, it must be cleaned, disinfected and sterilised before each change in patient.

Please ensure that only adequately validated device-specific and product-specific procedures are used for cleaning, disinfection and sterilisation, and that the validated parameters are observed in each cycle.

The cleaning and disinfection procedure must be recognised as effective (e.g., in the list of tested and recognised disinfectants and procedures issued by the Robert Koch Institute/DGHM) and already validated in all cases. If other procedures are used, their effectiveness must be confirmed during validation. CE-approved chemical disinfectants may also be used, provided they satisfy the requirements for material resistance (see Chapter 7.1 "Material resistance").

Please follow the instructions in Chapter 7.1 "Material resistance".

In addition, please follow the hygiene regulations of the hospital or medical practice.

6.4.1 Cleaning/disinfection

Clean and disinfect immediately after use. Ideally, a mechanical procedure (instrument dishwasher) is to be used.

- Prepare for cleaning and disinfection as described in Chapter 6.2 "Preparation".

Any suitable cleaning agent/disinfectant can be used.

Recommended procedure: Thermal disinfection

- Place the disassembled components of the PARI PEP System in the instrument disinfector.
- Select the 93 °C programme (treatment time 10 minutes).

The effectiveness of this procedure has been confirmed using a G7836 CD disinfector manufactured by Miele with neodisher®medizym cleaning agent manufactured by Chemische Fabrik Dr. Weigert, Hamburg.

or

Recommended procedure: Manual cleaning/disinfection

The effectiveness of this procedure has been proven using the Korsolex®-Endo-Cleaner cleaning agent in combination with the Korsolex®basic disinfectant made by BODE CHEMIE HAMBURG.

- The components must be placed in a 0.5% solution of Korsolex®-Endo-Cleaner and warm water at a temperature of approx. 50 °C, and cleaned thoroughly for 5 minutes.
- The components are then placed in a 4% disinfectant solution of Korsolex®basic for 15 minutes.
- Then the components of the PARI PEP System are rinsed thoroughly in hygienically tested water and left to dry completely (for at least 4 hours) on a clean, dry and absorbent surface.

Possible alternative: Chemical-thermal disinfection

- Place the disassembled components of the PARI PEP System in the instrument disinfector.
- Add a suitable cleaning agent/disinfectant and select the 60 °C programme.

Storage

- Between use, wrap the PARI PEP System in a clean, lint-free cloth and keep it in a dry, dust-free environment protected from possible sources of contamination.

6.4.2 Sterilisation

After cleaning/disinfection, wrap the dismantled PARI PEP System and the PEP tubing adapter if applicable in sterilisation packaging (disposable sterilisation packaging, e.g., foil/paper sterilisation bags). The sterilisation packaging must conform to DIN EN 868/ ISO 11607 and be suitable for steam sterilisation. Then, sterilise the PARI PEP System using the following sterilisation procedure:

Recommended procedure: Steam sterilisation

Sterilisation temperature: 121 °C (to be maintained for at least 20 min.) or
132 °C/134 °C (to be maintained for at least
5 min.)
max. 137 °C.

This procedure has been validated according to DIN EN 554/ISO 11134 as well as DGKH guidelines for the validation and routine monitoring of sterilisation processes using damp heat for medical devices.

Storage

Store the sterilised PARI PEP System in a dry, dust-free location away from sources of contamination.

7 Miscellaneous

7.1 Material resistance

Do not expose the PARI PEP S system to temperatures above 137 °C. Like all plastic parts, the PARI PEP System is susceptible to a certain degree of wear when it is used and reprocessed frequently. Over time, this can cause a change, which in turn impairs the therapeutic efficiency. Therefore, PARI GmbH recommends replacing the PARI PEP System, and the tubing adapter where applicable, at least once every two years when it is used regularly.

Note the following when selecting cleaning agents/disinfectants:

- Essentially, the aldehydic group of cleaning agents/disinfectants is the group that is suitable for cleaning and/or disinfecting this PARI PEP System.
- The use of other cleaning agents or disinfectants has not been tested with regard to the material resistance of this PARI PEP System.

7.2 Storage and transportation conditions

Do not store in damp rooms (e.g., bathrooms) or transport together with damp objects.

Protect from continuous direct sunlight during storage and transportation.

In general, the PARI PEP System should be stored as described in Chapter 6 "Hygienic re-use", "Storage" between applications. See Chapter 6.3 (at home) or Chapter 6.4 (hospital and doctor's practice).

8 Explanation of symbols

The following symbols can be found on the product and/or the packaging:

	Manufacturer
	The product satisfies the basic requirements as set out in Appendix I of Directive 93/42/EEC concerning medical devices.
	Please follow the instructions for use.
	Order no.
	Lot no.

Subject to technical changes.

For all product information please contact our Service Center:

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