
Zirela®

nebuliser handset


REF 678G8015

CE 0123

eFlow®
TECHNOLOGY

de Gebrauchsanweisung	fr Notice d'utilisation
en Instructions for use	it Istruzioni per l'uso
da Brugsanvisning	nl Gebruiksaanwijzing
es Istrucciones de uso	pt Manual de instruções
fi Käyttöohje	sv Bruksanvisning

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5 Cleaning and disinfection

The dismantled nebuliser handset (including the aerosol head) must be cleaned immediately after each use and disinfected at least once a day (additional information can be found in the instructions for use of the eFlow®/rapid nebuliser system).

6 Replacing the nebuliser handset

Each 28-day pack of Quinsair® contains a Zirela® nebuliser handset. Use the new nebuliser handset for each treatment cycle.

7 Package contents

One Zirela® nebuliser handset with one Zirela® aerosol head. There are no accessories for the Zirela® nebuliser handset.

8 Symbols

MD	Medical device
i	Follow the instructions for use
1	Single patient – multiple use
LOT	Lot number
UDI	Unique Device Identifier
REF	Item number
SN	Serial number
	Minimum and maximum ambient temperature
	Minimum and maximum humidity
	Minimum and maximum air pressure
CE 0123	CE marking: This product complies with the European requirements for medical devices (EU 2017/745).
	Do not dispose in domestic waste. This product must be disposed in accordance with the applicable (national, county or municipal) regulations.
	Caution
	Manufacturer
	Date of manufacture

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en Instructions for use

Zirela® nebuliser handset

- Read these instructions for use carefully and completely.
- Read the **Instructions for use of the eFlow®/rapid nebuliser system** completely before using for the first time. These contain further important instructions for the use of the Zirela® nebuliser handset, particularly regarding connection, operation, cleaning and disinfection.

Report serious incidents to the manufacturer and the competent authority.

1 Intended purpose

The Zirela® nebuliser handset with an eBase® Controller is intended to administer Quinsair® as an aerosol for inhalation.

To ensure correct dosing, the Zirela® nebuliser handset must only be used with Quinsair®.

2 Indication

The Zirela® nebuliser handset is indicated to administer Quinsair® to treat lung infections caused by *Pseudomonas aeruginosa* in patients with cystic fibrosis (CF).

We reserve the right to make technical changes.

Information as of: 2023-05

3 Contraindication

For the Zirela® nebuliser handset no contraindication is known to PARI Pharma GmbH.

Take note of contraindications of the medicine in the Quinsair® package leaflet.

4 Use of the nebuliser handset

The Zirela® nebuliser handset must only be used in accordance with its intended purpose. For reasons of hygiene and to avoid cross-infections, the Zirela® nebuliser handset must be used by one patient only.

- Before using for the first time, disassemble and clean the nebuliser handset including the aerosol head as described in the instructions for use of the eFlow®/rapid nebuliser system.

• Assemble the nebuliser handset as shown in the diagram:

- 1 Open the nebuliser handset.

- 2 Insert the aerosol head. **Do not touch the membrane (the curved metal surface in the middle of the aerosol head).**

- 3 Insert the inspiratory valve.

- 4 Close the nebuliser handset.

- 5 Attach the mouthpiece.

- 6 Connect the connection cord.

- 7 Fill with the medication; empty an ampoule completely.

- 8 Close the medication reservoir.

- 9 Sit in an upright position and relax.

This makes inhalation easier and improves uptake of the medication in the airways.

Hold the nebuliser handset horizontally all the time while performing inhalation. If the handset is tilted, this may result in the incorrect dosage being delivered.

- Hold the mouthpiece between your teeth and enclose it with your lips. The lips must not touch the blue expiratory valve.

The inhalation session should last about 5 minutes. Inhalation is ended as soon as no more aerosol is generated.

- Always check the medication reservoir to ensure that the entire dose of Quinsair® has been nebulised.

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