
Lamira®

Nebuliser System

eFlow
TECHNOLOGY

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en Instructions for use

da Brugsanvisning

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no Bruksanvisning

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sv Bruksanvisning


Kurzanleitung / Short guide

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For full details, please read and take note of
the instructions for use.

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1 ASSEMBLE

* TouchSpray made under licence from The Technology Partnership plc.

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
For full details, please read and take note of the instructions for use.

2 INHALE


Sit in an upright position.
Hold nebuliser handset horizontally.

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Warmes Wasser, Spülmittel

Clean immediately after inhalation

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
Disinfect

Boil thoroughly in distilled water for 5 minutes, or disinfect in the steam disinfecter for 6 minutes.

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Dry

Allow to air-dry completely.
Store in a clean, dry place.

de	Gebrauchsanweisung	3
en	Instructions for use	31
da	Brugsanvisning.....	59
es	Instrucciones de uso.....	87
fi	Käyttöohje	115
fr	Notice d'utilisation.....	143
it	Istruzioni per l'uso	173
nl	Gebruiksaanwijzing	201
no	Bruksanvisning.....	229
pt	Manual de instruções	257
sv	Bruksanvisning.....	287

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For Lamira® Type 678G2040

Table of contents

1	IMPORTANT INFORMATION	33
	Structure of the warnings	33
	Handling the nebuliser system.....	33
	Treating people who need assistance	34
2	PRODUCT DESCRIPTION.....	35
	Intended purpose.....	36
	Intended use.....	36
	Intended patient groups	36
	Indication(s)	36
	Contraindication(s)	36
	Lifetime of the medical device.....	36
3	GETTING STARTED	37
	Before using for the first time	37
	Preparing the controller	37
	Mains power	37
	Preparing the nebuliser handset.....	38
4	PERFORMING INHALATION	41
	Pausing during inhalation.....	42
	Screens during inhalation.....	42
	End of an inhalation session	43
5	CLEANING AND DISINFECTION	44
	Nebuliser handset.....	44
	Preparation	44
	Hygienic reprocessing at home	45
	Cleaning.....	45
	Disinfection	45
	Hygienic reprocessing in the hospital and doctor's practice	47
	Controller, connection cord and power adapter.....	47
	Disinfection	48
	Storage	48

6 TROUBLESHOOTING	49
7 DISPOSAL	52
8 SPARE PARTS AND ACCESSORIES	52
9 TECHNICAL DATA	52
General.....	52
Lamira nebuliser handset.....	53
Particle size distribution according to EN 13544-1	53
Classification according to IEC 60601-1 and IEC 60601-1-2	54
Materials used	54
Ambient conditions for operation	55
Ambient conditions for transport and storage.....	55
10 SYMBOLS	56
11 SERVICE AND CONTACT.....	57

1 IMPORTANT INFORMATION

Read these instructions for use and the instructions for use of the accessories carefully and completely before using them for the first time. Keep the instructions in a safe place in case you need to read them again later.

WARNING:

If the instructions for use are not followed, injuries and damage to the product cannot be ruled out.

Report serious incidents to the manufacturer and the competent authority.

Always seek the advice of a doctor before treating a disease.

Structure of the warnings

In these instructions for use, safety-related warnings are organised into hazard levels:

- The signal word **WARNING** is used to indicate dangers that may result in serious injuries or even death if precautionary measures are not taken.
- The signal word **CAUTION** is used to indicate dangers that may result in slight to moderate injuries or impair treatment if precautionary measures are not taken.
- The signal word **NOTE** is used to indicate general precautionary measures for handling the product that should be followed to avoid damaging the product.

Handling the nebuliser system

Check the nebuliser system before each use.

CAUTION:

Replace all broken, deformed or heavily discoloured parts. Damaged parts can impair the function of the nebuliser system and negatively affect the treatment.

WARNING:

Do not use the nebuliser system if the power adapter shows visible signs of damage. Otherwise there is a hazard from live parts (e.g., electric shock).

Therefore, follow the instructions below in order to operate the nebuliser system safely:

- Always unplug the power adapter from the socket to ensure complete disconnection from the mains.

WARNING:

Never disconnect the power adapter from the socket with wet hands. You might get an electric shock.

- Do not allow pets (e.g., rodents) to get close to the cables.

WARNING:

- Do not use the nebuliser system in areas where there is a risk of explosion or in the presence of gases that promote combustion (e.g., oxygen, nitrogen oxide, flammable anaesthetics).

- Keep the nebuliser system away from hot surfaces (e.g., hob).
- Do not use the system while operating motorised vehicles (risk of accident).

Treating people who need assistance


WARNING:

People who need assistance must be supervised constantly by an adult during the nebuliser treatment. This is the only way to ensure safe and effective treatment. Such persons often underestimate the hazards associated with the devices (e.g., risk of strangulation with the cable). This can lead to injuries.

WARNING:

The product contains small parts. Small parts can block the airways and pose a choking hazard. Therefore, make sure that they are always kept out of the reach of children.

2 PRODUCT DESCRIPTION

1 eBase® Controller

1a ON/OFF button

2 Lamira® nebuliser handset
(incl. aerosol head)

3 Lamira® Aerosol Head

4 Medication cap

5 Medication reservoir and

5a Aerosol chamber (preassembled)

6 Inspiratory valve

7 Mouthpiece with expiratory valve

8 Connection cord (connection between
the controller and the nebuliser
handset)

9 Power adapter

9a Country-specific adapter
(AU, EU, UK, US)

- Carrying case

- Batteries

Check that all components are contained in the package. If the package does not contain all the components, contact your service partner (see "SERVICE AND CONTACT" on page 57).

Intended purpose

The Lamira® Nebuliser System is intended to administer ARIKAYCE® liposomal 590 mg nebuliser dispersion as an aerosol for inhalation.

Intended use

The nebuliser system must only be used in accordance with its intended purpose.

WARNING:

For reasons of hygiene and to avoid cross-infections, the Lamira® Nebuliser System must be used by one patient only.

Only original PARI accessories must be used to operate the system.

The device is maintenance-free.

CAUTION:

Using accessories and spare parts other than those supplied by the manufacturer may cause increased electromagnetic emissions or reduced electromagnetic immunity, leading to a malfunction.

PARI Pharma GmbH is not responsible for damage to the system if it is used improperly or contrary to the intended use.

The device must only be serviced by PARI Pharma GmbH or a service location expressly authorised by PARI Pharma GmbH.

Intended patient groups

The nebuliser system is only suitable for use by patients who can breathe unassisted and are conscious.

The nebuliser system is intended for use by patients for whom ARIKAYCE® has been prescribed.

Indication(s)

The Lamira® Nebuliser System is intended to administer ARIKAYCE® for the treatment of non-tuberculous mycobacterial (NTM) lung infections in adults. For the indication(s) of ARIKAYCE®, please refer to the package leaflet of the medication.

Contraindication(s)

No contraindications are known for the Lamira® Nebuliser System. For contraindications of ARIKAYCE®, please refer to the package leaflet of the medication.

Lifetime of the medical device

- The expected lifetime of the eBase® Controller is 3 years.
- Each pack of ARIKAYCE® also contains a Lamira® nebuliser handset and 4 Lamira® aerosol heads. The nebuliser handset is designed for a treatment cycle of 28 days. Replace the aerosol head every week, **every time you start a new pack of seven vials.**

3 GETTING STARTED

Before using for the first time

- Wash your hands thoroughly every time before using the nebuliser system.
- To ensure that the nebuliser handset functions properly, ARIKAYCE must be brought to room temperature. Therefore, take ARIKAYCE out of the refrigerator 45 minutes before its application.
- In order to prepare and use ARIKAYCE, follow the instructions in the package leaflet for the medication.

WARNING:

Clean and disinfect the nebuliser handset and aerosol head before using them for the first time (CLEANING AND DISINFECTION, Page 44).

Preparing the controller


The controller can be operated with the delivered batteries or with commercially available rechargeable batteries.


Info:

Rechargeable batteries cannot be charged in the controller!

- Open the battery compartment by lifting the clip.

- Insert the (rechargeable) batteries in to the compartment so that they match the symbols +/-.

- Close the battery compartment.

Mains power


Make sure that the power adapter is undamaged before every use.

WARNING:


If the power adapter is damaged or faulty, do not use it.

The eBase® Controller was developed for mobile use. Therefore, insert (rechargeable) batteries before use and leave them in the controller during mains operation. This guarantees smooth functionality. This does not deplete the energy reserves of the (rechargeable) batteries.

- Replace the country-specific adapter for the power adapter as required. To do this, press on the clip and push the country-specific adapter upwards. When inserting the adapter (9a), make sure that it engages in the power adapter housing.

- Connect the power adapter to the power cord socket on the controller.

- Plug the power adapter into an electrical wall socket.

Preparing the nebuliser handset



CAUTION:

Check the nebuliser handset parts before every use. Replace all broken, deformed or heavily discoloured parts. Damaged parts and incorrectly assembled nebuliser handsets can impair the function of the nebuliser system and negatively affect the treatment.

- Insert the aerosol head (3) into the openings on the medication reservoir (5) (opened aerosol chamber (5a)). It must snap noticeably into place.


CAUTION:

Do not touch the membrane (the curved metal surface in the middle of the aerosol head) to avoid damage.


- Insert the inspiratory valve (6) in to the aerosol chamber.

When doing so, ensure that both valve flaps are seated correctly and face the aerosol chamber after insertion. They must lie perfectly flush.


- Close the aerosol chamber.


NOTE:

If the snap closure does not close, check whether the inspiratory valve has been inserted correctly. It must lie flush against the aerosol chamber.


- Fasten the mouthpiece (7) to the aerosol chamber.


- Attach the connection cord (8) to the nebuliser handset.


- Attach the connection cord (8) to the controller.


- Add the contents of a vial of ARIKAYCE to the medication reservoir.

- Close the medication reservoir.

- Check that all parts are connected tightly and that the medication reservoir is sealed.


4 PERFORMING INHALATION

- Sit in an upright position and relax. This makes inhaling easier and improves uptake of the medication in the airways.
- Hold the nebuliser handset in your hand.

CAUTION:


Hold the nebuliser handset horizontally all the time while performing inhalation. If the handset is tilted, this may result in the incorrect dosage being delivered.

- Press the ON/OFF button on the controller to start nebulisation.

A green LED lights up and an audible signal (1 beep) is emitted to indicate readiness for operation.

- Hold the mouthpiece between your teeth and enclose it with your lips. The lips must not touch the blue expiratory valve.

- Breathe in and out slowly and deeply through the mouthpiece. The mouthpiece should remain in your mouth while you are breathing out as well. Do not breathe through your nose. You can use a nose clip if necessary. The inhalation session should last approximately 14 minutes, but it can take up to 20 minutes.

Info:

It is normal for some fluid to collect in the aerosol chamber during inhalation. This is not a malfunction.

Info:

It is normal for aerosol to escape through the expiratory valve on the mouthpiece while you exhale. This is not a malfunction.

Pausing during inhalation

The inhalation session can be interrupted briefly if necessary by using Pause mode:

- Press the ON/OFF button for approx. 1 second.
Pause mode is activated when the LED light flashes green.
- Press the ON/OFF button for approx. 1 second again to resume the inhalation session.

Info:
Pause mode can only be activated after an operating time of a few seconds.

In the case of longer pauses, switch the controller off by pressing the ON/OFF button for at least 2 seconds.

Provided there is enough liquid in the medication reservoir, nebulisation can be restarted by pressing the ON/OFF button.


Screens during inhalation

The following indicators appear on the display during inhalation:


Start

During inhalation


Nebulising with battery power


Nebulising with mains power

Pause mode activated


Pause mode when operating with batteries

Pause mode when operating with mains power

End of an inhalation session

Inhalation is ended as soon as no more aerosol is generated. An audible signal (2 beeps) is emitted and the following indicator appears in the display of the controller:

The controller switches itself off automatically. If the controller switches itself off considerably sooner or later than usual, follow the instructions (see "TROUBLE-SHOOTING" on page 49).

Always check to make sure that the entire dose of ARIKAYCE has been nebulised.

- Open the medication cap and check the contents of the medication reservoir.
- If more than a few drops of ARIKAYCE are still in the medication reservoir, continue the inhalation session by putting the medication cap back in its place and then pressing the ON/OFF button.

Info:

When the system is used regularly, the time needed to complete nebulisation may slowly increase. This is not a malfunction.

- Operation with the power adapter:
When the inhalation session has ended, disconnect the power adapter from the socket.
- Carry out the hygiene procedures (see "CLEANING AND DISINFECTION" on page 44) immediately after the inhalation session has ended.

5 CLEANING AND DISINFECTION

Nebuliser handset

WARNING:

The nebuliser handset and the aerosol head must be cleaned and disinfected immediately after each use.

Preparation

Dismantle the nebuliser handset into its individual parts:


- Pull the connection cord out of the nebuliser handset.
- Remove the mouthpiece from the nebuliser handset.
- Open the medication reservoir.

Info:

The medication reservoir and aerosol chamber are already preassembled and must not be detached from one another.


- Open the aerosol chamber (5a).

- Remove the inspiratory valve from the aerosol chamber.

- Remove the aerosol head (3).

To do this, gently press the lateral engaging hooks together.

CAUTION:

Do not touch the membrane (the curved metal surface in the middle of the aerosol head) to avoid damage.

Hygienic reprocessing at home

WARNING:


Clean the nebuliser handset and the aerosol head immediately after each application.

NOTE:

Placing all the parts of the nebuliser handset and the aerosol head in warm washing water straight after use prevents medication residue from accumulating. This helps you to prevent a potential lengthening of the nebulisation times due to partial closure of the pores in the membrane.

Cleaning

- Rinse all nebuliser handset parts and the front and back of the aerosol head (for 10 seconds each) under warm running tap water (water of drinking water quality, temperature up to 37 °C).
- Place all nebuliser handset parts and the aerosol head in warm tap water with a little washing up liquid for 5 minutes. If the parts are heavily soiled, move the nebuliser handset parts gently back and forth in the washing water.

- Rinse all nebuliser handset parts and the front and back of the aerosol head thoroughly under warm running tap water to remove the dishwashing soap. Check each part, and if soiling is still visible on any part of the nebuliser handset or the aerosol head, place the part in washing water for another 5 minutes.

CAUTION:

Do not use any brushes or other tools to clean the aerosol head. The possibility of product components being damaged cannot be ruled out.

- Shake all parts to remove most of the water.

Disinfection

WARNING:

Disinfect the nebuliser handset and aerosol head immediately after cleaning them (the nebuliser handset cannot be disinfected effectively unless it has been cleaned beforehand).

NOTE:

Do not use a microwave to disinfect the aerosol head. The possibility of product components being damaged cannot be ruled out.

Boiling in distilled water


CAUTION:

Plastic may melt if it comes into contact with the hot bottom of the pot. Therefore, make sure there is enough water in the pot. Otherwise, the parts may be damaged.

- Place all parts in boiling water for at least 5 minutes. Use a clean pot and distilled water for this.

- Take the parts out of the pot.
- Shake all parts to remove most of the water.
- Place the individual parts on a dry, clean surface and allow them to dry completely.

CAUTION:

Moisture promotes the growth of bacteria. Please ensure adequate drying after every disinfection.

CAUTION:

Check the individual parts of the nebuliser handset after every disinfection. Replace all broken, deformed or heavily discoloured parts.

Possible alternative: Using a commercially available thermal disinfecter for baby bottles

- To ensure effective disinfection, use a thermal disinfecter with an operating time of at least 6 minutes. Information on how to carry out the disinfection, and for how long, and about the amount of water required is included in the instructions for use of the thermal disinfecter you use.

CAUTION:

Inadequate disinfection promotes the growth of bacteria and thus increases the risk of infection. If disinfection is not carried out for the full duration of the time specified for the respective thermal disinfecter, it is ineffective. Therefore, do not switch the device off before the disinfection cycle is complete. Also make sure that the thermal disinfecter is clean, and check it regularly to confirm that it functions properly.

- Take the parts out of the thermal disinfecter immediately after disinfection is complete. Place them on a clean, dry surface and allow them to dry.

CAUTION:

Check the individual parts of the nebuliser handset after every disinfection. Replace all broken, deformed or heavily discoloured parts.

Hygienic reprocessing in the hospital and doctor's practice

Use a washer-disinfector (WD) according to the requirements set out in the ISO 15883 series of standards.

The effectiveness of this procedure has been validated using the alkaline-enzymatic cleaning agent neodisher® MediClean forte (Chemische Fabrik Dr. Weigert, Hamburg) as the cleaning agent in conjunction with demineralised water as a neutraliser.

WARNING:

Carry out cleaning with disinfection immediately after use.

- Prepare for cleaning with disinfection (see "CLEANING AND DISINFECTION" on page 44).
- Place the nebuliser handset components in the washer-disinfector. Position the parts so that they are cleaned everywhere and no areas are left unwashed. The nebuliser handset is exclusively intended for use without change of patients and must not be used by several patients one after the other. Therefore, it is important to group nebuliser handset parts together when cleaning patients' devices together.
- Use demineralised water and a suitable cleaning agent, taking into account the instructions provided by the cleaning agent manufacturer.
- Select a cleaning and subsequent thermal disinfection programme that achieves an A0 value of min. 3000.

CAUTION:

Check the individual parts of the nebuliser handset after every disinfection. Replace all broken, deformed or heavily discoloured parts.

CAUTION:

Moisture promotes the growth of bacteria. Even when the washer-disinfector is equipped with a drying function, ensure that no residual moisture remains in the nebuliser handset.

- Remove any residual moisture by shaking out the parts, and place them on a dry, clean surface and allow them to air-dry completely.

Controller, connection cord and power adapter

Clean the housing surface of the controller and the cable as necessary with a clean, damp cloth.

NOTE:

Do not hold the controller under running water and do not use liquid cleaning agents. If liquid gets into the controller, the electronics may be damaged and the controller might not work properly anymore.

If liquid does get into the controller, contact your service partner.

Disinfection

Disinfect the controller, the connection cord and the power adapter immediately after cleaning. Use an alcohol-based, standard commercial disinfectant wipe for disinfecting. To ensure safe application, follow the instructions for use of the disinfectant wipe. Effectiveness has been tested using Bacillol Tissues and Clinell Wipes. Allow the controller, the connection cord and the power adapter to dry completely.

Storage

Between applications, particularly during longer breaks in treatment, keep the nebuliser handset in a dry, dust-free place.





Pack the controller, power adapter and connection cord in the case provided for this purpose.

NOTE:

Leaking batteries can damage the controller. Therefore, take the (rechargeable) batteries out of the controller if you do not intend to use the nebuliser system for a longer period.

6 TROUBLESHOOTING

You can correct malfunctions that occur during daily operation yourself. The following table will help you to identify faults more easily and correct the causes.

Screen (flashes alternately)	Possible cause	Solution
	Cable connection to the controller or on the nebuliser handset interrupted.	Check cable connections.
	No medication added.	Fill with medication.
	Inhalation interrupted.	Resume inhalation by restarting the device (press ON/OFF button).
	The maximum operating time of 20 minutes per inhalation session has been exceeded.	Resume inhalation by restarting the device (press ON/OFF button).
	Battery charge level less than 50%.	Have new (rechargeable) batteries ready.
	Battery is flat.	Insert new batteries or connect the power adapter.

Fault	Possible cause/Solution
The controller cannot be activated (no beep, no orange or green light signal).	<ul style="list-style-type: none"> - To switch the device on, press and hold the ON/OFF button for approx. 2 seconds. - Is the power plug connected to the power supply and the device? - Are the (rechargeable) batteries inserted correctly? - Check the charge level of the (rechargeable) batteries.
After switching the controller on, no aerosol is generated, and/or the controller switches itself off again after a few seconds.	<ul style="list-style-type: none"> - Has the medication reservoir been filled? - Check the connection between the controller and the nebuliser handset.
The controller switches itself off even though there is still medication in the medication reservoir.	<ul style="list-style-type: none"> - Check the charge level of the (rechargeable) batteries. - The nebuliser handset was not being held horizontally. Hold the nebuliser handset horizontally. - The controller switches itself off automatically after 20 minutes. Press the ON/OFF button again to continue the inhalation session and inhale the prescribed dosing of medication.
The device does not switch off automatically even though no more aerosol is being generated and there is no more medication in the medication reservoir.	<p>The inhalation session has now been completed successfully and the device can be switched off by pressing the ON/OFF button.</p>

Fault	Possible cause/Solution
The nebulisation times have become noticeably longer.	<p>When the system is used regularly, the time needed to complete nebulisation may slowly increase. This is not a malfunction. However, if you notice that the nebulisation times have become significantly longer, this may be due to one of the following reasons:</p> <ul style="list-style-type: none"> - The medication cap was not closed correctly. - The aerosol head was not cleaned immediately after the last application, or it is damaged. - You are using the aerosol head for longer or more often than is recommended in section Lifetime of the medical device.
Aerosol escapes continuously in large volumes from the opening slots in the medication reservoir.	<ul style="list-style-type: none"> • Check whether the nebuliser handset has been assembled correctly. When doing so, ensure that both valve flaps of the inspiratory valve are seated correctly. When the valve is inserted, they must lie flat and face towards the aerosol chamber.
During inhalation, a small amount of fluid escapes through the mouthpiece.	<p>A small amount of fluid collects in the aerosol chamber during inhalation, and when the nebuliser handset is tilted it can escape through the mouthpiece.</p> <ul style="list-style-type: none"> • Hold the nebuliser handset horizontally.
A small amount of fluid drips from the bottom of the expiratory valve during the inhalation session.	<p>Some drops of fluid collect on the expiratory valve on the mouthpiece during treatment and can run down the mouthpiece.</p> <ul style="list-style-type: none"> • Keep a cloth to hand.
The medication cap cannot be taken off the medication reservoir after inhalation treatment.	<p>The negative pressure in the medication reservoir is too great.</p> <ul style="list-style-type: none"> • Open the aerosol chamber and carefully pull the aerosol head to enable air to flow in. <p>The medication cap can now be removed easily.</p>

If the malfunction persists, contact your service partner.

7 DISPOSAL

Live parts (such as the controller, cable, power adapter and (rechargeable) batteries) as well as the aerosol head must not be disposed of with normal domestic waste.

The device components must be disposed of in accordance with local regulations.

8 SPARE PARTS AND ACCESSORIES

Description	Item number
Nebuliser handset (incl. aerosol head)	678G8250
Power adapter (100 – 240 V~, 50-60 Hz)	078B7114
Connection cord	178G6009

9 TECHNICAL DATA

General

Electrical connection	PARI power adapter REF 078B7114: Input: 100 V – 240 V~, 50-60 Hz Output: 12 V ---
Batteries	4 x 1.5 V Mignon (AA LR6/Alkaline)
Rechargeable batteries	4 x 1.2 V Mignon (AA rechargeable, at least 2100 mAh)
Power consumption	Approx. 2 W
Weight of the controller (incl. batteries)	approx. 300 g
Sound pressure level (at a distance of 1 m)	35 dB +/- 1 dB
Dimensions of the controller housing	H 4.0 cm, Ø 11.6 cm

Technical electromagnetic compatibility data is available in table form upon request from PARI Pharma GmbH or on the internet at the following link:
<https://www.pari.com/fileadmin/Electromagnetic-compatibility-4.pdf>


Lamira nebuliser handset

Nebuliser type	Membrane nebuliser
Weight of the nebuliser handset	approx. 55 g
Dimensions of the nebuliser handset (W x H x D)	5.0 cm x 6.7 cm x 14.5 cm
Min. / Max. fill volume	8.4 ml
MMAD ^{a)} (Mass Median Aerodynamic Diameter)	4.4 µm
Particle spectrum ^{a)}	4.3 µm to 4.8 µm
Respirable fraction ^{a)} < 5 µm	55.7% +/- 3.6%
Aerosol output ^{b)}	4.46 ml
Aerosol output rate ^{b)}	0.35 ml/min

Aerosol data according to EN 13544-1; Alternative test method in accordance with Annex CC.
Nebulised medium: ARIKAYCE 8.4 ml.

- a) Measurement with Next Generation Pharmaceutical Impactor (NGI) at 23 °C and 50% relative humidity. Inspiratory flow: 15 l/min.
- b) Measurement with breathing simulator at 23 °C and 50% relative humidity. 500 ml breathing volume, breathing frequency 15 cycles/minute, sinusoidal breathing pattern, inhalation/exhalation ratio 1:1 (with adults, may vary in children).

Particle size distribution according to EN 13544-1

Classification according to IEC 60601-1 and IEC 60601-1-2

Type of electric shock protection (power adapter)	Protection class II
Degree of protection from electric shock of applied part (nebuliser handset)	Type BF
Degree of protection against ingress of water as per IEC 60529 (IP rating)	IP 21
Degree of protection when used in the presence of flammable mixtures of anaesthetics with air, with oxygen or with nitrous oxide	No protection
Operating mode	Continuous operation
Environment	- Professional healthcare facility environments - Domestic healthcare settings

Materials used

This product mainly consists of the following materials:

- Controller:
Thermoplastic elastomers, acrylonitrile butadiene styrene copolymer, polyethylene
- Nebuliser handset (incl. aerosol head):
Polypropylene, thermoplastic elastomers, polyoxymethylene, silicone rubber, stainless steel
- Mouthpiece:
Polypropylene, silicone rubber
Lamira® contains no parts made from natural rubber (latex).

Ambient conditions for operation

- Ambient temperature:
+5 °C to +40 °C
- Relative humidity of ambient air:
15% to 93% (non-condensing)
- Air pressure: 700 hPa to 1060 hPa

Within professional healthcare facility environments, the device must only be used in in-patient wards and intensive care units. The device must never be used in areas exposed to elevated magnetic or electrical radiation (e.g., in the vicinity of an MRI scanner). The device must not be used in close proximity to HF surgical equipment.

Users can operate the device both at home and when out and about. Its use on trains and aeroplanes is restricted to passenger areas only. If using the device in a car, it must be operated with battery power.

CAUTION:

When using portable wireless communication devices, you must ensure that a distance of at least 30 cm is maintained between all parts of the nebuliser system and these devices (including their accessories, e.g., antenna cables or external antennas). Otherwise, the performance of your nebuliser system may deteriorate.

CAUTION:

The controller must not be placed directly beside or stacked with other devices for operation. If operation near or stacked with other devices is unavoidable, the controller must be monitored continuously to ensure proper operation.

Ambient conditions for transport and storage

- Temperature: -25 °C to +70 °C
- Relative humidity of ambient air:
0% to 93% (non-condensing)
- Air pressure: 500 hPa to 1060 hPa

These values for transport are valid for the unpackaged product / the unpackaged aerosol head. The blister packing for the aerosol head will only function properly under the transport conditions described on the system packaging.

NOTE:

Moisture from condensation can impair the controller's functional capabilities. Extreme temperature fluctuations should therefore be avoided. Wait until the controller has adjusted to room temperature before you start an inhalation session.

If the device is not used for a prolonged period, the (rechargeable) batteries should be removed first.

Info:

Batteries vary widely with regard to life cycle and performance depending on their quality. The operating time of about 90 minutes can only be achieved with the batteries (4x) included with the delivery.

10 SYMBOLS

MD	Medical device
	Follow the instructions for use
	Single patient – multiple use
	Direct current
	Alternating current
LOT	Lot number
UDI	Unique Device Identifier
REF	Item number
SN	Serial number
	Protection class of the applied part: Type BF
	Protection class II equipment
IP21	The device is protected against solid foreign objects of 12.5 mm Ø and greater and against vertically falling water drops.
	Minimum and maximum ambient temperature
	Minimum and maximum humidity
	Minimum and maximum air pressure
CE 0123	CE marking: This product complies with the European requirements for medical devices (EU 2017/745).
	Do not dispose in domestic waste. This product must be disposed in accordance with the applicable (national, county or municipal) regulations.
	Caution

Manufacturer

Date of manufacture

11 SERVICE AND CONTACT

In case of technical problems or questions about the device, please contact your service partner. The contact details of your responsible service partner as well as a digital version of your Lamira® instructions for use can be found at www.pari.com/lamira.

We reserve the right to make technical changes.

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