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
## PARI O-PEP

- Model: PARI O-PEP (Type 018)

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
- Model: PARI O-PEP (Type 018)

## **Read the instructions for use**

Read these instructions carefully before using the product. Follow all instructions and safety directions. Keep the instructions in a safe place.

## **Validity of instructions for use**

PARI O-PEP (Type 018)

## **Contact**

Email: info@pari.de  
Tel.: +49 (0)8151-279 220 (international)  
+49 (0)8151-279 279 (German)

## **Copyright**

©2025 PARI GmbH Spezialisten für effektive Inhalation  
Moosstrasse 3, 82319 Starnberg, GERMANY. Technical and design changes and printing errors reserved. The same applies to illustrations.

## **Disclaimer**

These instructions for use describe the components of PARI products and optional accessories. For this reason, these instructions for use also describe and illustrate features not present in your PARI product because they are, for instance, country-specific and/or optional. When using the systems, products and functions, the applicable country-specific regulations must be observed.

## **Trade marks**

Registered trade marks of PARI GmbH Spezialisten für effektive Inhalation in Germany and/or other countries:

PARI®

# TABLE OF CONTENTS

<b>IMPORTANT INFORMATION</b> .....	<b>4</b>
Intended purpose.....	4
Indication.....	4
Contraindication.....	4
Labelling.....	5
Safety and warning instructions.....	5
<b>PRODUCT DESCRIPTION</b> .....	<b>7</b>
Components.....	7
Overview and designations.....	7
Description of function.....	8
Material information.....	8
Operating life .....	8
<b>USE</b> .....	<b>9</b>
Assembling the PARI O-PEP.....	9
Perform treatment .....	10
<b>REPROCESSING</b> .....	<b>12</b>
Reprocessing cycles.....	13
Preparation for reprocessing.....	13
Reprocessing the PARI O-PEP.....	14
Reprocessing the case and neck strap.....	16
Further validated processes for reprocessing.....	17
<b>FURTHER INFORMATION</b> .....	<b>17</b>

# 1 IMPORTANT INFORMATION

## 1.1 Intended purpose

The PARI O-PEP (**O**scillating **P**ositive **E**xpiratory **P**ressure) is used for the mobilisation of secretion in cases of acute and chronic diseases of the lower airways.

The PARI O-PEP can be used:

- by children aged from 5 to 8 under the supervision and direction of a competent person
- by children aged 8 and older and adults after informed instruction

Individuals in the following groups must be supervised by a person who is responsible for their safety:

- Infants and children
- Individuals with limited capabilities (e.g., physical, mental, sensory)

If the patient is not able to use this product safely on their own, then the treatment must be carried out by the responsible person.

This PARI product is designed for patients who are able to breathe by themselves and are conscious.

PARI O-PEP device is provided non-sterile. It is re-usable with defined cleaning and disinfection procedures.

This PARI product can be used in a home environment, as well as in professional health institutions. When used in a home environment, this PARI product is intended for single-patient use only (no patient change). In a professional environment, the device can be used with different patients as long as the corresponding hygiene reprocessing measures are complied with.

The frequency and duration of use is determined by professional medical staff<sup>1</sup> according to the individual needs of the patient.

## 1.2 Indication

Diseases of the lower airways which are associated with increased mucus build up.

## 1.3 Contraindication

The PARI O-PEP must not be used by people who are suffering untreated pneumothorax or massive haemoptysis.

---

1) Professional medical staff: Doctors, pharmacists, and physiotherapists.

## 1.4 Labelling

The following symbols can be found on the product and/or the packaging:

	Medical device
	Unique Device Identifier (UDI)
	Legal manufacturer
	Date of manufacture
	Serial number
	Item no.
	Production batch number, lot number
	This product conforms to the EU Medical Device Regulation 2017/745.
	Consult instructions for use
	O-PEP
	Steel ball

## 1.5 Safety and warning instructions

**The present instructions for use contain important information, safety instructions and precautionary measures. The user must follow these in order to guarantee safe operation of this PARI product.**

This PARI product must be used only as described in these instructions for use.

### Labelling and classification of warning instructions

In these instructions for use, safety-critical warnings are categorised according to the following hazard levels:

**DANGER**

DANGER indicates a hazardous situation which will lead to very severe injuries or death if it is not avoided.

## **WARNING**

WARNING indicates a hazardous situation which can lead to very severe injuries or death if it is not avoided.

## **CAUTION**

CAUTION indicates a hazardous situation which can lead to mild or moderate injuries if it is not avoided.

## **NOTE**

NOTE indicates a hazardous situation which can lead to material damage if it is not avoided.

## **General**

The PARI O-PEP must not be used unless instruction has been given by a healthcare professional<sup>1</sup>. The PARI O-PEP cannot be used together with a mask. Masks have exhalation valves through which the air would escape during exhalation and would therefore not enter the PARI O-PEP. Respiratory therapy is therefore not possible if a mask is worn.

**If your health condition does not improve or it even worsens as a result of the treatment, seek professional medical advice.**

## **Hazard due to small parts which can be swallowed**

### **DANGER**

#### **Risk of suffocation due to blocked airways**

The product contains small parts. Small parts can block the airways and lead to a choking hazard.

- Keep all components of the product out of the reach of babies and infants at all times.

## **Reporting serious incidents**


Report serious adverse incidents to the manufacturer and to the competent authority.

## 2 PRODUCT DESCRIPTION

### 2.1 Components

The product components included in the scope of delivery are country-specific, and the contents may deviate from the components described in this instruction for use. For the scope of delivery, therefore please refer to the packaging.

### 2.2 Overview and designations

1	Lower part with mouthpiece
2	Funnel
3	Steel ball
4	Upper part
5	Neck strap
6	Case

## 2.3 Description of function

When the patient repeatedly exhales into the device against a mobile steel ball, a vibration is created which is transferred to the lungs. This vibration loosens the mucus in the lower respiratory tract. During use, the respiratory tract is kept open, and the removal of the loosened mucus is improved.

If the PARI O-PEP is kept horizontal, the oscillation frequency during use is approx. 15 Hz. This frequency can be modified by tilting the device upwards (higher frequency) or downwards (lower frequency). This allows selection of the appropriate frequency for each patient.

## 2.4 Material information

The individual product components are made from the following materials:

Product component	Material
Lower part with mouthpiece, funnel, upper part	POM (polyoxymethylene)
Steel ball	Stainless steel
Neck strap	Polyester
Case	Polyester

## 2.5 Operating life

Product component	Max. service life	Processing limits
PARI O-PEP	1 year	300 disinfection cycles, 100 sterilisation cycles

The product's service life depends on the reprocessing and inhalation cycles. It may be reduced through repeated reprocessing or multiple daily use. The product components are designed for the cycles specified. The individual components should be replaced no later than the maximum service life or the specified reprocessing limits.

## 3 USE

All the steps described below must be carried out properly.

Do not use product components unless they have been thoroughly cleaned and dried. Wash your hands thoroughly before every use. You must perform cleaning and disinfection before using the device for the first time.

### 3.1 Assembling the PARI O-PEP

#### WARNING

##### **Danger of injury through strangulation**

For individuals who are not able to perform the therapy session without assistance or cannot appreciate the hazards, the risk of injury is greater, e.g. strangulation by the neck strap. Such individuals include, for example, babies, children, and people with limited capabilities.

- Ensure that for these individuals a person responsible for their safety either supervises or implements the application.

#### WARNING

##### **Choking hazard from breathing in small parts**

The product is not completely closed. As a result, there is a danger that small parts can get inside it, and may be breathed in during inhalation.

- Therefore, before every use it is imperative to make sure there are no foreign bodies inside the product.


#### CAUTION

##### **Risk of impaired treatment/injury**


Damaged components and/or an incorrectly assembled nebuliser may impair the functioning of the product and therefore the treatment.

- Check all product components and the accessories before each use.
- Replace any broken, deformed or seriously discoloured parts.
- Follow the assembly instructions in these instructions for use.

1. Attach the lanyard to the lower part.

1. Insert the funnel in the lower part.
2. Place the ball in the funnel.
3. Screw the upper part onto the lower part.

### 3.2 Perform treatment

**All the safety instructions and warnings in these instructions for use must have been read and understood before any treatment is carried out.**

The standard treatment using the PARI O-PEP is described below. It is possible that your doctor or therapist recommends a usage which differs. In order to achieve the desired therapeutic effect using the PARI O-PEP, always perform the treatment in the way recommended by your doctor or therapist.

#### CAUTION

##### **Danger of overexertion**

Overexertion during use can cause dizziness. This can happen more quickly with children, because they have a smaller lung capacity than adults.

- If dizziness occurs, discontinue the treatment immediately.
- Before the next treatment, speak to your doctor or therapist.

1. Use the neck strap to hang the PARI O-PEP around your neck.
2. Sit in an upright position and relax.
3. Take the PARI O-PEP in your hand, and hold it so that the mouthpiece is aligned horizontally.
4. Hold the mouthpiece between your teeth and enclose it with your lips.

---

 Ensure that the holes in the upper part remain clear, and are not covered by your hand.

---

5. Breathe in as slowly and deeply as possible through your nose.
6. Hold your breath for approx. 1 to 2 seconds.
7. Breathe out slowly and deeply through your mouth into the PARI O-PEP (as slowly and long as possible). For this, tension your cheek muscles, and keep your cheeks taut.
8. Repeat breathing in and out for 10 to 15 breaths. End the treatment when the exertion becomes too great. Rapid, strained breathing renders the treatment ineffective.

---

***Instructions for optimum treatment:***

 While inhaling and exhaling, place your free hand flat on your chest. During correct use, you can sense a clear, throbbing vibration. The strength of the vibration changes with the alignment of the PARI O-PEP. Find the alignment at which the throbbing is the strongest (tilt the PARI O-PEP somewhat further upwards or downwards).

---

## 4 REPROCESSING

### WARNING

#### Risk of infection

When products are reused, there is an increased risk of bacterial transmission and infections, especially through cross-contamination, bacterial growth, and residual organisms or moisture

- Be absolutely sure to perform cleaning and disinfection, even before using the device for the first time.
- To avoid residual organisms and cross-contamination, clean, disinfect and, if necessary, sterilise all products as per the reprocessing instructions and the specified reprocessing cycles.
- To avoid bacterial growth caused by residual moisture, it is essential that you ensure complete drying before sterilisation, and also after reprocessing.
- Always use clean reprocessing materials, and ensure that the product is cleaned fully. Contamination which remains affects further reprocessing.

### CAUTION

#### Danger of infection for patients at risk

For patients at risk, airway infections represent an increased risk of deterioration of general health, because these patients are particularly vulnerable to residual germs. Patients at risk include cystic fibrosis patients, patients with immunosuppression or immunodeficiency, and vulnerable patient groups.

- If you are a patient at risk, disinfect the individual parts once per day during the therapy.
- If you are unsure whether you are a patient at risk, consult with specialist medical personnel before use.

---

*Always use drinking water for reprocessing in a home environment.*

 *In professional health care institutions, drinking water is sufficient for precleaning. For all other reprocessing steps and processes, use deionised water with a low microbiological load (at least drinking water quality).*

---

Please observe the instructions for use for the chemicals used.

Wash your hands thoroughly before every reprocessing.

Medical aids for reprocessing must be clean.

## 4.1 Reprocessing cycles

Perform the reprocessing steps as per the specified reprocessing cycles.

### Reprocessing cycles in a home environment

PARI O-PEP (incl. funnel and steel ball)	– Clean immediately after every use – Disinfect once a week
Neck strap, case	Cleaning visible dirt

### Reprocessing cycles in professional healthcare institution

#### *SINGLE PATIENT USE*

PARI O-PEP (incl. funnel and steel ball)	– Clean immediately after every use – Disinfect once a week
Neck strap, case	Cleaning visible dirt

#### *BEFORE A CHANGE OF PATIENTS*

PARI O-PEP (incl. funnel and steel ball)	– Cleaning – Disinfection – Sterilisation
Neck strap	– Replacement
Case	– Cleaning visible dirt – Wipe disinfection

## 4.2 Preparation for reprocessing

### Procedure

<b>PARI O-PEP</b>	Dismantle the PARI O-PEP into its individual parts:
	1. Unscrew the upper part from the lower part.
	2. Remove the steel ball from the funnel. Make sure that the steel ball does not roll away.
	3. Remove the funnel from the lower part.
	4. Remove the neck strap from the lower part.


### 4.3 Reprocessing the PARI O-PEP

The following products can be cleaned, disinfected and sterilised according to the procedure described below:

– PARI O-PEP

The neck strap and the case [see: Reprocessing the case and neck strap, page 16] must be treated separately.

#### Procedure

<b>Step 1: Pre-cleaning</b>	All individual parts must be precleaned immediately after use. Avoid residues on the device immediately after use. – Tap water at a temperature $\leq 40\text{ }^{\circ}\text{C}$ 1. Immediately after use, rinse all individual parts (except case and neck strap) under flowing tap water, in order to remove visible contamination.	
<b>Home environment</b>	<b>Step 2: Cleaning</b> – Tap water at a temperature $\leq 40\text{ }^{\circ}\text{C}$ – Standard commercial washing-up liquid – Container having sufficient capacity 1. Unless otherwise specified by the manufacturer of the washing-up liquid, add about 1 teaspoon of washing-up liquid to 3 l warm drinking water. 2. Place all the parts in the washing-up water. Application time: 5 minutes 3. Occasionally move the parts back and forth. 4. In case of visible soiling, use a medium-soft brush (e.g. a toothbrush) which is used exclusively for this purpose. 5. Remove the individual parts, and rinse them (except case and neck strap) under flowing tap water, in order to remove soap residues.	
<b>Home environment</b>	<b>Step 3: Disinfection</b> To ensure a sufficient level of disinfection, you must clean the product beforehand. <b>Thermal disinfection in boiling water</b> – Clean cooking pot – Drinking water (recommendation: use deionised water instead of this) 1. Place all the individual parts in water at a rolling boil for at least 5 minutes. 2. Ensure that there is sufficient water in the saucepan.	

## Procedure

<b>In professional health care institutions</b>	<b>Step 2: Cleaning and disinfection</b>	<ul style="list-style-type: none"> <li>– Cleaning and disinfection device as per DIN EN ISO 15883-1 and 15883-2</li> <li>– Suitable push-in trolley with instrument holder</li> </ul> <hr/> <p><i>As a rule, automatic disinfection forms part of a cleaning and disinfection device programme which also includes automatic cleaning, e.g.:</i></p> <ul style="list-style-type: none"> <li>– Pre-rinse with cold tap water for 2 minutes.</li> <li> – Wash at 55 °C for 5 minutes, using 0.5 % alkaline-enzyme cleaner.</li> <li>– Perform an intermediate rinse for 5 minutes, using cold demineralised water.</li> <li>– Thermal disinfection at 90 °C in contamination-free water (&lt;10 cfu/ml) for 5 minutes (A0&gt;3000).</li> </ul> <hr/> <ol style="list-style-type: none"> <li>1. Position the individual parts in the instrument holder in such a way that they are thoroughly cleaned all over.</li> <li>2. Start the programme.</li> </ol>
	<b>Drying</b>	<p>Ensure that there no residual moisture remains in the components after the cleaning or disinfection. Dry the device as follows:</p> <ol style="list-style-type: none"> <li>1. Shake the water out of all of the parts.</li> <li>2. Place all parts on a dry, clean and absorbent surface and allow them to dry completely.</li> </ol>
<b>In professional healthcare institutions (before patient change)</b>	<b>Step 3: Sterilisation</b>	<p>To ensure sufficient sterilisation, you must clean and disinfect the product beforehand.</p> <p>For sterilisation, use a process compliant with the ISO 17665 series of standards.</p> <ul style="list-style-type: none"> <li>– Temperature: 134 °C</li> <li>– Holding time: minimum 3 minutes to maximum 5 minutes</li> </ul> <ol style="list-style-type: none"> <li>1. Pack all the disassembled parts in a sterile barrier system according to ISO 11607-1 (e.g. foil-paper packaging).</li> <li>2. Sterilise the product as per ISO 17665, using a steam steriliser at the following parameters: Fractionated pre-vacuum at 134 °C for between 3 and 5 minutes.</li> </ol>
<b>Visual inspection</b>	<ol style="list-style-type: none"> <li>1. Inspect all product components after each cleaning, disinfection or, where applicable, sterilisation. Replace defective, deformed or seriously discoloured parts.</li> </ol>	

## Procedure

### Storage

Store this product as described below:

- Dry and dust-free, e.g. in a clean, lint-free cloth (e.g. a tea towel)
  - Protected against contamination, where necessary (e.g. using optional sterile packaging)
1. Push the steel ball into the small side compartment.
  2. Re-assemble the PARI O-PEP (except steel ball) [see: Assembling the PARI O-PEP, page 9].
  3. Place the PARI O-PEP, together with the neck strap, in the case.

## 4.4 Reprocessing the case and neck strap

*The neck strap cannot be disinfected or sterilised.*

## Procedure

### Step 1: Cleaning

1. Use a clean, damp cloth to wipe the outer surface of the case and the neck strap.

### Step 2: Disinfection (before patient change)

For disinfection, use a standard, alcohol-based disinfectant (e.g. isopropanol). When applying and dosing the disinfectant, it is essential to follow the instructions for use for the product.

1. If there is visible contamination, clean the case before disinfection.
2. Moisten a cloth with the disinfectant.
3. Wipe the outer surface of the case thoroughly with the cloth.
4. Let the disinfectant dry completely.

## 4.5 Further validated processes for reprocessing

The instructions provided were validated by PARI and were found to be suitable for preparing your medical device for its reuse.

Further validated procedures for reprocessing:

[https://www.pari.com/fileadmin/041D0624\\_Professional\\_healthcare\\_institution\\_Validated\\_Reprocessing\\_Methods.pdf](https://www.pari.com/fileadmin/041D0624_Professional_healthcare_institution_Validated_Reprocessing_Methods.pdf)

Ensure that the reprocessing actually performed by your personnel on your equipment with the chemicals used achieves the desired results. For this, validation and routine monitoring of the process are usually required. In particular, if you have to deviate from our validated process, ensure that the reprocessing method selected by you is appropriately effective, and that potential adverse effects are assessed.

## 5 FURTHER INFORMATION

### 5.1 Disposal

All product components may be disposed of with normal domestic waste. The country-specific disposal regulations must be observed.

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PARI GmbH, Spezialisten für effektive Inhalation