

Quality - Performance - Innovation



PRONEB[®] *Ultra II* COMPRESSOR NEBULIZER SYSTEM

The PARI PRONEB[®] Ultra II delivers the quality and performance you've grown to trust from PARI in a smaller and lighter package. Featuring:



PARI LC[®] Reusable Nebulizers

Used in more clinical trials by leading Pharmaceutical companies than all other brands combined.

PARI LC[®] SPRINT #23F35

PARI LC[®] PLUS #22F81

Small, Light, and Compact

Specialists in effective inhalation



PRONEB® Ultra II

COMPRESSOR NEBULIZER SYSTEM

The PARI PRONEB® Ultra II Compressor comes in three configurations and choice of LC® Reusable Nebulizer:



PRONEB® Ultra II Deluxe

Standard compressor configuration plus an instructional video DVD and carrying case.

PARI LC® SPRINT VERSION #86F35-LCSD

PARI LC® PLUS VERSION #86F81-LCD



PRONEB® Ultra II Pediatric

Featuring our green pediatric compressor and PARI Pals™ sticker set. Also includes the instructional DVD, carrying case, backup nebulizer, and Bubbles the Fish™ II Pediatric Mask.

PARI LC® SPRINT VERSION #86F35-P

PARI LC® PLUS VERSION #86F83-P



PRONEB® Ultra II Standard

Grey compressor, PARI LC® Reusable Nebulizer, and backup nebulizer.

PARI LC® SPRINT VERSION #86F35-LCS

PARI LC® PLUS VERSION #86F83-LC+



PRONEB® ULTRA II CARRYING CASE #44F2212

ADULT PVC MASK #14F20-PVC

BUBBLES THE FISH™ II PEDIATRIC AEROSOL MASK #44F7248

PARI PALS™ STICKER SET #ML-44-734

PARI BABY™ CONVERSION KIT
SIZE 1 (1-12 MTHS) #44F2401
SIZE 2 (1-3 YRS) #44F2402
SIZE 3 (3+ YRS) #44F3301



PRONEB® Ultra II Compressor with PARI LC® Reusable Nebulizers



PARI LC® SPRINT REUSABLE NEBULIZER



PARI LC® PLUS REUSABLE NEBULIZER



PARI LC® STAR REUSABLE NEBULIZER

	TOR	MMD	%RF	ATT
PARI LC® SPRINT	500 mg/min	3.5 µm	68%	~5 min
PARI LC® PLUS	440 mg/min	3.8 µm	65%	6-7 min
PARI LC® STAR	400 mg/min	3.1 µm	78%	6.5-8 min

TOR= TOTAL OUTPUT RATE [mg/min]

MMD= MASS MEDIAN DIAMETER¹ [µm]

%RF= PERCENTAGE RESPIRABLE FRACTION [0.1~5µm]

ATT= AVERAGE TREATMENT TIME [min]

¹Measured with Malvern MasterSizer X at 50% relative humidity, 0.9% NaCl solution, inspiratory flow 20 liters/minute, continuous nebulization, 23° C, fill volume 2.5 ml.



Long-Lasting Air Filter

Easy to install and remove
AIR FILTER #41E485IPI2

PRONEB® ULTRA II Compressor Technical Specifications:

Maximum Pressure of Compressor	34 PSI (2.3 bar)
Maximum Liter Flow of Compressor	8.5 LPM
Nominal Operating Pressure ²	18psi (1.2 bar)
Nominal Operating Air Flow ²	4.5 LPM
Weight	3 lbs (1300 gms)
Dimensions (Length, Width, Height)	6.5 x 6.5 x 3.75 in
Warranty	5 years

²When used with a PARI LC® Reusable Nebulizer

Specialists in effective inhalation



Call PARI at 1.800.FAST.NEB (327.8632) to Order or Contact Your Wholesaler