



PARI PEP FAMILY

POSITIVE **EXPIRATORY** PRESSURE



PARI PEP FAMILY



For patients diagnosed with Cystic Fibrosis or lung disease with secretory problems, PARI PEP devices help efficiently loosen mucous and clear the airways.

PEP THERAPY THREE WAYS

PARI O-PEP
Oscillating Therapy



PARI PEP™ S
Combined Therapy



PARI PEP™ S
Conventional Therapy



Configurations

Part

PARI PEP™ S

018F63

PARI O-PEP

018F65

Accessories

Part

PARI LC® Sprint

023F35

PARI LC PLUS®

022F81

Noseclip

041F3500

Pressure Monitor

044B4612

Not for use with PARI O-PEP

*PARI PEP systems approved for ages 5 and above

MEET THE PEP FAMILY

PARI O-PEP OSCILLATING THERAPY



Configuration includes the PARI O-PEP oscillating device, lanyard, carrying case, and instructions for use.

The PARI O-PEP is an Oscillating Positive Expiratory Pressure device in which patients exhale into the device against a movable ball, causing a vibration that is transmitted to the lungs. The O-PEP produces a range of oscillation frequencies between 6 and 20 Hz and can be adjusted by changing the incline of the device.

PARI PEP™ S - COMBINATION PEP/NEBULIZER THERAPY OR TRADITIONAL PEP THERAPY

Option 1:

Combined PEP and PARI LC® Reusable Nebulizer Therapy



Option 2:

Conventional PEP therapy



Configuration includes the PARI PEP™ S resistance ring and mouthpiece, inspiratory valve for conventional therapy, and instructions for use.

Designed for use with the PARI LC PLUS® or LC® Sprint nebulizers, combination therapy provides concurrent PEP therapy with clinically proven aerosol delivery. PEP™ S offers eight resistance settings covering a large expiratory flow range that offers pressure of 10-20cm H₂O.



PARI Respiratory Equipment, Inc.
Tel: 1-800-FAST-NEB (327-8632)
product.info@pari.com
www.pari.com | Midlothian, VA 23112