

PARI Vios® delivers consistent particle size with targeted aerosol delivery to the lungs in just 5-6 minutes.¹ Available in both standard and pediatric configurations, PARI Vios includes environmentally friendly packaging, instructions, and 100% recycled cotton tote to help reduce our carbon footprint.

## **ADULT AND PEDIATRIC OPTIONS**



Standard Configuration (Includes Back Up Nebulizer: not shown)



Pediatric configuration also includes Bubbles the Fish™ II mask, carry tote, and activity set

Configurations		
	Part #	HCPCS #
Standard, LC® Sprint	310F35-LCS	E0570, A7005
Standard, LC PLUS®	310F83-LC+	E0570, A7005
Pediatric, LC <sup>®</sup> Sprint	310F35-P	E0570, A7005, A7015
Pediatric, LC PLUS®	310F83-P	E0570, A7005, A7015

Technical Specifications		
Operating Pressure <sup>2</sup>	18 PSI (1.2 bar)	
Liter Flow <sup>2</sup>	4.5 LPM	
Weight	3 lbs	
Dimensions (L, W, H)	6.5 x 6.5 x 3.75 in	
Compressor Warranty	5 years	
1 Massured with Malvern MasterSizer V at E00/ relative humidity, 0.00/ NaCl solution		

- <sup>1</sup> Measured with Malvern MasterSizer X at 50% relative humidity, 0.9% NaCl solution, inspiratory flow 20 liters/minute, continuous nebulization, 23° C, fill volume 2.5 ml.
- <sup>2</sup> When used with a PARI reusable nebulizer.

## PATIENTS BENEFIT FROM:



Fast 5-6 minute treatments help improve patients' quality of life



Breath-enhanced nebulizer design saves money by wasting less medication



Consistent particle size with targeted medication delivery to the lungs



PARI LC Sprint 023F35



PARI LC PLUS 022F81

## **Accessories**

Air Filter



Bubbles the Fish™ II Pediatric Mask 044F7248



PARI Adult Aerosol Mask 044F7252

